# USER MANUAL BEDIENUNGSANLEITUNG

# **WINNER**

English	5
Deutsch	47

# **English**

### © 2008 Handicare

### All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorisation of Handicare.

The information provided is based on general data concerning the constructions known at the time of the publication of this manual. Handicare executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. Handicare cannot be held liable for possible damage resulting from specifications of the product deviating from the standard configuration.

The available information has been prepared with all possible diligence, but Handicare cannot be held liable for possible errors in the information or the consequences thereof. Handicare accepts no liability for loss resulting from work executed by third parties.

Names, trade names, etc. used by Handicare may not, as per the legislation concerning the protection of trade names, be considered as being available.

2008-01

### **Table of contents**

Intro	duc	on	. 8
	Deliv	ery	. 8
	This	nanual	. 8
	Indic	ations and pictograms	. 9
	Picto	grams	11
	Avai	able documentation	11
	Serv	ce and technical support	11
	Iden	fication platefication plate	12
	Iden	fication of the product	12
	War	anty stipulations concerning the scooter	13
		ity stipulations concerning the Product	
	Use	scooters and the environment	15
	Use	according to purpose1	15
1	Gen	ral safety regulations and instructions	
	1.1	Decals and instructions on the scooter	16
	1.2	Technical specifications	16
	1.3	Modifications	16
	1.4	Safety	16
2	Gen	ral description	19
3	Maiı	components	20
A	A		
4	•	ation	
4	<b>Upe</b> 4.1	Driving with the Winner	21
4	4.1	Oriving with the Winner	21 21
4	4.1	Oriving with the Winner	21 21 22
4	4.1	Oriving with the Winner	21 21 22 22
4	4.1	Driving with the Winner	21 21 22 22 25
4	4.1 4.2	Oriving with the Winner	21 22 22 25 25
4	4.1 4.2	Driving with the Winner	21 22 22 25 25 25
4	4.1 4.2	Driving with the Winner  4.1.1 Inspection before use  Deration  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection  Driving  4.3.1 Driving forward and in reverse	21 22 22 25 25 25 26
4	4.1 4.2	Driving with the Winner  4.1.1 Inspection before use  Deration  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection  Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners	21 21 22 22 25 25 26 27
4	4.1 4.2	Driving with the Winner  4.1.1 Inspection before use  Degration  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection  Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners  4.3.3 Slopes	21 22 22 25 25 25 26 27
4	4.1 4.2	Driving with the Winner  4.1.1 Inspection before use Degration  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection  Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners  4.3.3 Slopes  4.3.4 Obstacles	21 21 22 25 25 25 26 27 27 28
4	4.1 4.2	Driving with the Winner  4.1.1 Inspection before use Degration  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection  Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners  4.3.3 Slopes  4.3.4 Obstacles  4.3.5 Overload	21 21 22 25 25 25 26 27 27 28 29
4	4.1 4.2	Driving with the Winner	21 21 22 25 25 25 27 27 28 29
4	4.1 4.2	Driving with the Winner  4.1.1 Inspection before use Degration  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection  Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners  4.3.3 Slopes  4.3.4 Obstacles  4.3.5 Overload  4.3.6 Switching off  4.3.7 Parking	21 22 22 25 25 25 27 27 28 29 29
4	4.1 4.2	Driving with the Winner  4.1.1 Inspection before use.  Degration.  4.2.1 Dashboard.  4.2.2 Switching on.  4.2.3 Charging connection.  Driving.  4.3.1 Driving forward and in reverse.  4.3.2 Corners.  4.3.3 Slopes.  4.3.4 Obstacles.  4.3.5 Overload.  4.3.6 Switching off.  4.3.7 Parking.  4.3.8 Charging the batteries.	21 22 22 25 25 25 27 27 28 29 29 29
4	4.1 4.2 4.3	Driving with the Winner	21 21 22 25 25 25 27 27 28 29 29 30
4	4.1 4.2 4.3	Driving with the Winner  4.1.1 Inspection before use Degration  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners  4.3.3 Slopes  4.3.4 Obstacles  4.3.5 Overload  4.3.6 Switching off  4.3.7 Parking  4.3.8 Charging the batteries  4.3.9 Battery charger  Transferring in and out	21 22 22 25 25 25 27 27 29 29 29 30 31
4	4.1 4.2 4.3	Driving with the Winner.  4.1.1 Inspection before use.  Operation	21 22 22 25 25 25 27 27 29 29 31 31 32
4	4.1 4.2 4.3 4.4 4.5	Driving with the Winner	21 22 22 25 25 27 28 29 29 29 31 31 32
4	4.1 4.2 4.3	Driving with the Winner  4.1.1 Inspection before use Departion  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners  4.3.3 Slopes  4.3.4 Obstacles  4.3.5 Overload  4.3.6 Switching off  4.3.7 Parking  4.3.8 Charging the batteries  4.3.9 Battery charger  Transferring in and out  Pushing the scooter  4.5.1 Neutral position of the scooter  Reducing the Winner	21 22 25 25 25 27 27 29 29 31 32 32 33
4	4.1 4.2 4.3 4.4 4.5	Driving with the Winner  4.1.1 Inspection before use Departion  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners  4.3.3 Slopes  4.3.4 Obstacles  4.3.5 Overload  4.3.6 Switching off  4.3.7 Parking  4.3.8 Charging the batteries  4.3.9 Battery charger  Transferring in and out  Pushing the scooter  4.5.1 Neutral position of the scooter  Reducing the Winner  4.6.1 Removing the seat	21 22 22 25 25 26 27 28 29 29 31 32 33 33 33
4	4.1 4.2 4.3 4.4 4.5 4.6	Driving with the Winner.  4.1.1 Inspection before use.  Deparation.  4.2.1 Dashboard.  4.2.2 Switching on.  4.2.3 Charging connection.  Driving  4.3.1 Driving forward and in reverse.  4.3.2 Corners.  4.3.3 Slopes.  4.3.4 Obstacles.  4.3.5 Overload.  4.3.6 Switching off.  4.3.7 Parking.  4.3.8 Charging the batteries.  4.3.9 Battery charger.  Transferring in and out.  Pushing the scooter.  4.5.1 Neutral position of the scooter.  Reducing the Winner.  4.6.1 Removing the seat.  4.6.2 Folding down the steering column.	21 22 22 25 25 26 27 28 29 29 31 32 33 33 33
4	4.1 4.2 4.3 4.4 4.5	Driving with the Winner  4.1.1 Inspection before use Departion  4.2.1 Dashboard  4.2.2 Switching on  4.2.3 Charging connection Driving  4.3.1 Driving forward and in reverse  4.3.2 Corners  4.3.3 Slopes  4.3.4 Obstacles  4.3.5 Overload  4.3.6 Switching off  4.3.7 Parking  4.3.8 Charging the batteries  4.3.9 Battery charger  Transferring in and out  Pushing the scooter  4.5.1 Neutral position of the scooter  Reducing the Winner  4.6.1 Removing the seat	21 22 22 25 25 27 29 29 31 32 33 33 33 34

5	Adj	ustments options	35
	5.1	Seat adjustments	35
		5.1.1 Adjusting the seat height	
		5.1.2 Adjusting the seat depth	
		5.1.3 Backrest adjustment	
		5.1.4 Lower-backrest adjustment	
		5.1.5 Adjusting the armrest	36
		5.1.6 Adjusting the headrest (option)	
	5.2	Adjusting the steering column	37
6	Mai	ntenance	38
	6.1	Maintenance table	38
		6.1.1 Batteries	38
		6.1.2 Replacing the batteries	39
		6.1.3 Cleaning the batteries	40
	6.2	Tyres	40
		6.2.1 Inflating the tyres	40
	6.3	Cleaning	41
7	Tro	ubleshooting	42
	7.1	Troubleshooting table	42
	7.2	Automatic fuse	43
	7.3	EMC malfunction	43
8	Tec	chnical specifications	45
	8.1		
		8.1.1 Winner	45
	8.2	Approval	46
	8.3	CE-Declaration	46
	8.4	Authorised service and technical support	46

### Introduction

Welcome to a growing community of users who use a scooter from the Handicare product range.

Handicare stands for reliability and innovative technology, which results in easily operated quality products.

### **Delivery**


Your dealer will deliver the Winner scooter ready for use.

### This manual

With this manual you can safely use and maintain (clean) the product. In case of doubt always contact your dealer.

In the documentation the words "left", "right", "front" and "back" are used to indicate a certain part of the product. This is viewed from the position of the user.

### Indications and pictograms

The product has the following decals and labels.

- A Charging connection
- B Fuse for lighting
- C Automatic fuse
- D Freewheel
- E Tyre pressure wheels
- F Adjusting the steering column

### A. Charging connection

For charging of the batteries, see 'Charging'.


### B. Fuse for the lighting

Here you find the fuse for the lighting. See 'Product specifications sheet'.

### C. Automatic fuse


The scooter has a safety fuse feature. See 'automatic fuse' and 'faults'.

### D. Freewheel


- 1. Lever in the lowermost position:
  - The drive of the motor concerned is disconnected: the wheelchair can be pushed.
- 2. Lever in the uppermost position:

The drive of the relevant motor has been switched on: The wheelchair can be electrically driven.

### E. Tyre pressure wheels

For tyre pressure of the wheels, see 'product specification sheet'.

### F. Adjusting the steering column

Adjust the steering column to obtain a position that is most comfortable for you while driving. See 'Adjusting the steering column'.

The steering column may also be folded down to reduce the size of the scooter for transport purposes. See 'Folding down the steering column'.

### **Pictograms**

The following pictograms are used in this manual:

### Caution

Procedures that can result - if they are not executed properly - in damage to the product, the environment, or human injury.

### Attention!

Suggestions and advice that help tasks or procedures involved to be carried out more easily.

Consult the specified information source(s) first.

Pull the charge cord out of the charging connection of the scooter before doing any maintenance on the scooter.

### **Available documentation**

The following technical documentation is available for this scooter:

- User manual
- Service manual.

### Service and technical support

For information concerning specific settings, maintenance or repair work please contact your dealer. He is always prepared to help you.

Ensure you have at hand:


- Model
- Year of manufacture
- Identification number

This information is provided on the identification plate. See 'Identification of the product'.

### **Identification plate**

See Fig. 00000010.

On this plate (A) information on the product can be found. See: 'Identification of the product'.


### Identification of the product

See Fig. 010007010.

- A Model
- B Year of manufacture
- C Identification number
- D User area indoors or outdoors
- E Maximum load in kg

For the place of the identification plate, see 'main components'.

### Warranty stipulations concerning the scooter

In the following warranty and liability stipulations the terms and definitions as summed up hereafter must be explained as follows:

Product: The hand-operated or electric wheelchair or scooter

manufactured and delivered by Handicare.

Customer: The person who directly obtains a Product from Handicare or an

authorised representative.

Dealer: The person who delivers a Product obtained from Handicare to

customers or third parties.

User: The person who uses a Product manufactured by Handicare.

Irrespective of what is determined concerning warranty conditions applicable to the Product, in any case the following applies with regard to the warranty:

- 1. Except insofar as described otherwise Handicare guarantees the Product for its suitability for the purpose for which the Product is intended all of these points as described in this manual and for the quality of the material of which the Product is made and the manner in which the Product is manufactured.
- 2. Repairs or replacement of parts of the Product that may be necessary as a result of faults that are based on qualitatively faulty material or manufacturing errors will be executed free of charge, as long as such faults occurred within one (1) year after the date of delivery of the Product to the Customer. The parts to be replaced must be shipped postage-paid to Handicare. Installation or removal of these parts is for the Customer's account. Therefore, the following cases are not eligible for free repair or replacement as meant above:
  - The repair or replacement that is necessary in connection with faults that arise after one (1) year after the date of delivery of the Product to the Customer;
  - The repair or replacement that is required in connection with faults due to improper or careless use of the Product or that are based on the Product being used for another purpose than the one for which it is intended, in which regard it applies that if the Customer is a Dealer, this Dealer will indemnify Handicare against possible claims from Users or other third parties for faults based on an incorrect or careless use of the Product;
  - Parts that are subject to wear, and the need for repair or replacement of the
    parts are the actual consequence of normal wear. These can include, but are
    not limited to, tyres, shrouds, arm-pads, seatings etc.
- 3. Irrespective of that stipulated under 2, it applies as far as an electric Product is concerned that with regard to the battery that forms part of the Product warranty is only given in case of faults or non-functioning of the battery that are demonstrably the direct consequence of material or manufacturing errors. A fault or non-functioning of the battery as a result of normal wear is not covered by the warranty as meant in these warranty stipulations. Similarly not covered by the warranty are faults or non-functioning that are the consequence of improper or unprofessional use of the Product or the battery that is part of the Product, including the incorrect charging of the battery and the failure to carry out timely and good maintenance, in which context it also applies that if the Customer is a Dealer, this Dealer will indemnify Handicare against possible claims from Users or other third parties that are based on the above mentioned improper or unprofessional use of the Product

- or the battery that is part of the Product. Handicare does not guarantee batteries provided by dealers or purchased by the customer.
- 4. The warranty conditions as mentioned in above stipulations in any case become null and void if:
  - The guidelines of Handicare for the maintenance of the Product are not, or insufficiently, followed;
  - A necessary repair or replacement of parts is based on neglect, damage or abuse of the Product or a use of the Product for another purpose than the one for which it was intended;
  - Parts of the Product are replaced by parts of another origin than those used by Handicare and/or parts of the Product are replaced without the permission of Handicare.
- 5. The warranties as mentioned in stipulations 1 through 3 also become null and void if what is involved is re-used by a new user within the warranty period and that such re-use necessitates adaptations to the product and those adaptations are carried out without the instructions and/or at the order of Handicare.
- 6. To retain rights under the warranties delineated above, the Customer must, in case of damage or other calamities, as rapidly as possible contact Handicare and inform them as fully as possible. The possibility of taking recourse to the above-mentioned warranty conditions becomes null and void for the Customer in any case after 20 workdays after the claim situation or the calamity arises that was the reason for the recourse to the warranty.
- 7. The replacement of a part or the repair or the reconditioning of the Product within a current warranty period does not extend the warranty period.
- 8. Handicare gives no warranty on repair to or reconditioning of the Product carried out other than under order of and/or at the instructions of Handicare. If repairs and/or reconditioning are executed by or on behalf of a Customer, the Customer indemnifies Handicare with respect to the claims of third parties who result, in the broadest sense of the word, from such repairs or reconditioning.

### **Liability stipulations concerning the Product**

Irrespective of what is determined regarding liability in the general conditions applicable to the Product, with regard to liability in any case the following applies:

- taking into consideration the following stipulations, Handicare only accepts liability
  for loss due to death or physical injury that is the result of a defect in the Product for
  which Handicare is responsible and for damage to another object that is the private
  property of the user of the Product, as long as said loss is the direct result of a fault
  in the Product.
- 2. Handicare accepts no other or further liability than delineated under 1. In particular Handicare accepts no liability for consequential damage, in any form whatsoever.

### Used scooters and the environment

If your scooter has become superfluous or needs to be replaced, it can usually be taken back by your dealer. If this is not possible, please contact your local authorities for the possibilities of recycling or an environmentally friendly way of disposing of the used materials. For the production of the scooter several plastics and metals have been used. The scooter also contains electronic components that must go to electronic waste. Used batteries belong to chemical waste.

### Use according to purpose1

The scooter Winner was designed for:

- transport of persons with a weight of up to 159 kg
- use on paved roads, sidewalks, footways, and bike-ways.
- use in and around the house

Your supplier should give you good user's instructions before you can independently operate the product and participate in traffic.

You must be able to correct the consequences of actions while driving the Winner. Consider operating the Winner scooter under the monitoring of an experienced dealer during first operations.

The Winner scooter is not a motorised vehicle in the traffic regulations sense.

The Winner has a maximum speed of 15 km/h. (Winner Alpine 9 km/h)

# Please make sure you are well acquainted with the contents of this owner's manual before driving the winner!

If you use the scooter for other purposes than those for which it is intended, Handicare accepts no liability whatsoever for damage or injury resulting from such use other than that for which the scooter was developed and designed.

'Use according to purpose' as established in EN 292-1 is the use for which the technical product is suitable according to the statement of the manufacturer, including his instructions in the sales brochure. In case of doubt this is the use that follows from the construction, execution and function of the product. Use according to purpose also includes taking into consideration the instructions in the user manual.

### 1 General safety regulations and instructions

Handicare accepts no liability whatsoever for damage or injury caused by the failure to (strictly) adhere to the safety guidelines and instructions or else due to carelessness during the use and cleaning of the scooter and possible accessories. Depending on the specific operating circumstances or the accessories used, additional safety instructions may be required. Please contact your dealer immediately if you observe any potential danger during the use of the product.

The user of the scooter (see under 'use according to purpose') is at all times fully responsible for the fulfilment of the locally applicable safety regulations and guidelines.

### 1.1 Decals and instructions on the scooter

Signs, symbols and instructions placed on this scooter comprise part of its safety facilities. They must never be covered or removed and must remain present and clearly legible throughout the entire life span of the scooter.

Immediately replace or repair illegible or damaged indications, symbols and instructions. Contact your dealer with regard to this.

### 1.2 Technical specifications

The technical specifications must not be altered.

### 1.3 Modifications

Modifications to (parts of) this product are not permitted.

### 1.4 Safety

To avoid accidents and undesirable situations, it is of great importance to pay attention to the following safety instructions.

### Pay extra attention to driving on slopes:

- Never drive the Winner up a slope with an angle greater than 10°.
- On slopes always drive slowly and with great awareness.
- Never drive down a slope at full speed.
- Do not drive down slopes with loose gravel or a sandy surface, because one of the back wheels might slip.
- Do not turn on a slope.

- Never take a turn at full speed. Decrease speed when taking a turn.
- Emergency Brake: If a calamity occurs while driving, and you need to make a sudden stop, you must release the operating handle immediately.

- Be sure that no items of clothing hang loose. These could become caught between the wheels.
- Make sure your fingers do not get caught in the mechanism of the steering adjustment.

- Adjust your driving according to the circumstances:
  - Drive carefully on slippery roads, resulting from rain, ice, or snow!
  - Drive more slowly in crowded environments.
  - Do not use the scooter on uneven roads.
  - Avoid the Winner from coming into contact with sea water. Sea water is caustic and may damage the scooter.
  - Avoid the Winner from coming into contact with sand. Sand can affect the moving parts of the scooter, resulting in unnecessarily rapid wear.
- Never operate the scooter when you are under the influence of drugs, alcohol or medication that could affect your driving ability.
- You must have sufficient visual acuity to be able to drive safely in the scooter.
- It is imperative to turn on the lights in case of reduced visibility.
- The scooter has been equipped with brake lights that light up when you brake.
- Only use the indicator lights to indicate an intended change in driving direction.
- Only use the horn when it is necessary to warn pedestrians of other road users of dangerous circumstances.

- Never place metal parts on top of the batteries: this may cause a short circuit in the batteries, resulting in damage.
- Do not take along any passengers: your scooter was specifically designed to transport only yourself.
- Do not burden your scooter with any more load than can be placed in the basket.
- Never hang anything from the handlebars.
- Never use your scooter to pull a cart: the scooter was not designed for this. Anyway, it is prohibited to do so. You might seriously damage your scooter if you do.
- Do not use the foot plate as a step.
- Do not place your feet close to the front wheel and the front fork. Always keep your feet on the rubber mat, rather than against the front cover.

 The scooter is electronically operated. Certain parameters have been factory-set. These settings guarantee a comfortable and efficient use. It is not permitted to change these parameters.

- The standard version of your scooter was tested according to the strictest EMC requirements. Mobile telephones have no influence on the driving behaviour of the scooter.
   When using a mobile telephone in the vicinity of a scooter with special adaptations, you are advised to first switch off the scooter.
- Your scooter may influence electromagnetic fields such as alarm systems.
- If the electronics of the scooter is not properly screened off, this
  may affect sensitive electrical appliances, such as shop alarms
  and garage openers. The scooter has been tested for this.
  Should any problems of this nature occur, please report these to
  your supplier immediately.

- Take care in ultra-violet light: this may cause premature deterioration of materials such as rubber, plastic, and enamel.
- Watch out for protruding parts: these may damage your scooter or the environment.
- Keep your scooter away from open flames.

### 2 General description

The Winner is an electrically driven scooter for use outside, for longer distances.

The height of the seat is adjustable according to the leg length of the user. The seat can be turned, and the armrests folded aside, by the user when he/she wants to make a transfer.

The steering console is adjustable to enable a comfortable seating position.

When driving at higher speeds, one should be extra careful. Slow down, especially on sidewalks and in pedestrian areas.

The following should be kept in mind when using the scooter:

- A scooter is not a motorised vehicle in the traffic regulations sense. When participating in traffic, you have the same rights and duties as a pedestrian.
- A driver's license is not required. You do not have to pay any road tax. You do, however, need a motorbike third-party insurance.
- Using the Winner is not restricted to any age.
- At all times, you are responsible for keeping your scooter in perfect working order, so that driving the scooter is a safe act.
- The scooter is a technical product. Therefore, any repairs and maintenance should be carried out by qualified personnel. If any work needs to be done, you should contact your dealer.

### 3 Main components

See Fig. 000000020.

The following main components are, or may be, installed on the scooter:

A. Seat. This is where the user is seated.B. Backrest. To support the back of the user.

C. Head rest (option). To support the head.D. Armrest. To support the arms.

E. Steering column On this, the dashboard and all operating keys and handles

are located.


F. Lighting. When driving in the dark the lights must be on.G. Indicators. To indicate a change in direction while driving.

H. Detachable basket. You can carry your groceries in this.

I. Front wheel. This is also the driving wheel.

J. Front bumper. Protects the front wheel in case of a collision.K. Back wheel. Both back wheels are driven and suspended.

Depending on the application and the design, parts can be adjusted for optimum seating comfort, and parts can be removed for transportation. Refer to the sections concerned in this owner's manual.

### 4 Operation

### 4.1 Driving with the Winner

To be able to drive, this owner's manual describes the following parts.

- Inspection before use.
- Dashboard.
- Driving.
- Transferring in and out.
- Pushing.
- Reducing.

### 4.1.1 Inspection before use

Before one starts driving the scooter, it is important to check the following points of the scooter:

- make sure the seat is properly interlocked.
- make sure all clamp connections of the various adjustments are securely tightened, to prevent loosening by vibration.
- make sure all lights and indicators are working properly, both at the front and the back. switching on the switches on the dashboard, to make sure the lights and indicators are working properly.
- make sure the hand brake on the front wheel is released.
   Pull the handle towards you as far as possible to release the hand brake. This will cause the catch to be released. Now you can let go of the handle; the hand brake will be released.
- make sure the tyres are sufficiently pumped up. See 'Product specifications'.
  - Poorly pumped up tyres have a negative effect on the driving characteristics.
  - Poorly pumped up tyres cause a higher power consumption, depleting the batteries more quickly.
  - Poorly pumped up tyres causes unnecessary wear to the tyres.
- make sure the batteries are charged sufficiently: this is indicated by the green zone on the battery-condition meter. See 'Dashboard'.

In winter, batteries have a lower capacity. In case of light frost, the capacity drops to about 75% and at a temperature of less than -5 degrees, it will drop to about 50% of the normal capacity. This will decrease the range of action.

- make sure the scooter is in the driving position. See 'Neutral position of the scooter'.
- make sure the brakes are functioning correctly. See 'driving with the scooter'.

Should the scooter, despite all safety precautions, start behaving in an unpredictable manner, let go of the throttle handle and, after having come to a standstill, take the contact plug from the contact.


If everything is in order, the scooter can be driven.

### 4.2 Operation

See Fig. 000000030.


The scooter is equipped with the following operating commodities:

- A Dashboard with all operating switches
- B Operating handles forward / reverse
- C Fuse for the lighting. See 'Lighting fuse'
- D Charger connection. See 'Charging the batteries'.
- E Operating handle for steering-column adjustment
- F Front brake with parking lock

### 4.2.1 Dashboard

See Fig. 000000040.

The dashboard is provided with the latest technology and offers reliable and useful functions for the operation of your scooter.

The dashboard features the following organs / switches:

- A Indicator left and right\*
- B Battery-condition meter
- C Speed control
- D Lighting
- E Plug contact
- F Horn\*
- G Alarm lights
- \* These keys have been mounted both on the right and on the left of the dashboard.

### A. Indicators left / right

By operating the switch (on the left or on the right), the indicator concerned will start flashing, to indicate that you intend to change your driving direction:

- move the switch to the left if you wish to turn to the left;
- move the switch to the right if you wish to turn to the right.

As soon as you have completed your left or right turn, move the switch to its neutral position to switch off the indicator.

### B. Battery-condition meter

The battery-condition meter gives a general indication of the condition of the batteries. It is quite normal for the pointer to go down during acceleration, because the battery voltage is reduced when power is requested. This voltage reduction is no indication for the capacity of the batteries. When the batteries are fully charged, the pointer will be at its maximum position. However, this will also be the case if the batteries are not fully charged, owing to the battery's characteristics.

The most reliable battery indication is shown while driving on a flat surface.

The battery-condition meter has three areas in which the pointer may rest after the scooter is switched on.

When the battery is sufficiently charged, the pointer is in the green area. As the battery is depleted, the pointer will move to the yellow, and, subsequently, into the red area.

### The fields have the following meaning:

- Red: The batteries are depleted and have to be recharged as soon as possible.
- Yellow: The batteries are quite depleted.
   The scooter may be used, but the batteries should very soon be recharged.
- Green: The batteries are fully charged.

### C. Speed control

With this knob, the desired maximum speed is set. Turning clockwise will increase the speed (the maximum has been set by the manufacturer to 15 km/h).

Turning counterclockwise will reduce the maximum speed.

Select the position of this button beforehand, in accordance with the driving conditions, such as driving in cramped spaces or in a room with many people.

### D. Lighting

Pressing this switch will turn on the lighting:

The switch has the following positions:

0: lighting is off

1: lighting is on.

### E. Horn

The horn buttons are used to give a warning signal in dangerous situations. The horn will sound for as long as you keep the button pressed.

# F. Contact

To switch on the scooter, the contact plug has to be inserted into the contact. All of the Winner's functions, with the exception of the alarm lights and the lighting, only operate if the contact plug has been inserted into the contact.

The alarm lights and the lighting will always work, even if the contact plug has not been put into the contact.

### G. Alarm lights

Pressing this switch will turn on the alarm lights.

You should use the alarm lights if you fear you will not be seen by other traffic, or if you are stopped because of a malfunction.


Pressing this switch a second time will switch off the alarm lights again.

### 4.2.2 Switching on

See Fig. 00000050.

The scooter should be switched on as follows:

 Put the contact plug into the contact, as far as it will go (A).


### 4.2.3 Charging connection

See Fig. 000000060.

In the steering column, underneath the dashboard, there is the charger connection (A) to which the battery charger cable can be connected.

All electronics should be switched off during the charging of the batteries.

### 4.3 Driving

Before driving with the scooter, it should be adjusted optimally to your personal needs. Then, after you have carried out all check-ups, you can seat yourself in the scooter. See 'Transferring in and out'.

Now you may start your ride.

As a scooter driver, you are very vulnerable in traffic. Keep in mind that other traffic participants may not always notice you. Keep to the current traffic regulations and the Road Traffic Act.

Avoid out-of-the-way routes, so that help can be on its way quickly when needed.

# 4.3.1 Driving forward and in reverse

See Fig. 000000070.

Switch on the scooter. See 'Switching on'.

The operating handles are coloured:

Green: forward driving Yellow: reverse driving

### Forward driving

- Driving forward using thumb control: Slowly push forward the operating handle (A) with your thumb. The further you push the handle forward, the faster you will drive.
- Driving forward using hand control: Slowly pull back the operating handle (D) by hand. The further you pull the handle towards you, the faster you will drive.

- Reverse driving:
   Release the operating handle. The Winner will come to a standstill.
- Driving in reverse using thumb control:
   When the scooter has comet to a complete standstill, slowly push the operating handle (C) forward with your thumb. The further you push the handle forward, the faster you will drive in reverse.
- Driving in reverse using hand control:
   Slowly pull back the operating handle (B) by hand. The further you pull the handle towards you, the faster you will drive in reverse.

- Make sure there is nothing behind the scooter when you reverse.
- The maximum speed in reverse is equal to half the maximum speed forward.

### Braking and stopping while driving forward or reversing

If you release the operating handle slowly, the Winner will brake and come to a standstill.


While braking, the brake lights will go on.

Under normal circumstances, you have to release the operating handle slowly when you want to stop.

As soon as the scooter is stopped, the automatic parking brake is applied.

### Remarks:

You can inspect the brakes by suddenly releasing the operating handle.

- When driving at a high speed, you need to be extra careful, especially on sidewalks and in pedestrian areas. In that case, it is advisable to set a lower maximum speed using the speed controller.
- Do not cross a street until you are able to handle the scooter and its controls well.

### 4.3.2 Corners

### Moving around corners in forward direction

Turning left or right

By turning the handlebars to the left or to the right, the scooter will change direction and start moving to the left or to the right accordingly.

- You should always use your indicators when making a turn.
- Do not release the operating handle when making a turn, because this will bring the scooter to a standstill.
- Always take turns at a safe, often lower, speed.

### Reversing while making a turn

Turning left or right

By turning the handlebars to the left or to the right, the scooter will change direction and start moving in the opposite direction of the direction you steer.

- Steering to the right: the scooter will reverse while making a left turn.
- Steering to the left: the scooter will reverse while making a right turn.

### 4.3.3 Slopes

Maximum slope which can safely be driven with maximum weight according to ISO 7176-2 is 10° at 159 kg.


- Slopes with an angle greater than 10° may not be taken.
- Avoid sudden and jolting movements, such as sudden braking, to avoid the scooter becoming unstable.
- Do not brake too hard when moving around corners.
- Avoid changing directions or turning around on a slope.
- Do not turn on slopes steeper than 10°.

### Driving up a slope

Driving up a slope is done in the same way as driving in forward direction.

Lean slightly forward with the upper part of your body, while driving up a slope. This will provide the scooter with more stability.

- Drive up a slope at half speed. Drive at a constant speed.
- If you notice a strong decrease in speed when driving on a slope, choose a less steep route, to avoid the scooter from getting overheated.
- If you drive on a slope too long, this may result in overheating.
- If you want to drive up a steep slope and the speed of your scooter decreases drastically, you should find an easier route.

### Remarks:

Whenever possible, choose a road that is not so steep.

### Driving down a slope


If you need to drive down a slope, drive at the lowest possible speed and with the utmost control.

- Turn the speed control fully to the left (lowest speed) before driving down a slope.
- Move the upper part of your body slightly backward and always brake slowly. If you brake fast, you may cause the scooter to tip over forward.

### 4.3.4 Obstacles

See Fig. 000000090.

### **Driving up obstacles**

- Drive straight up to the kerb, almost touching it with your front wheel.
- Push the operating handle and drive onto the sidewalk without changing direction.
- As soon as the front wheel is on the sidewalk, you need to keep up your speed for the back wheels to be able to move up onto the sidewalk as well.
  - If you cannot drive up onto the sidewalk, find a lower part where your can.
  - Any obstacle should be taken in forward direction, and with a certain run-up speed.

The scooter has the power to take obstacles of 10 cm.

Practice taking slightly elevated obstacles. Once you are practised at these, gradually increase the elevation of the obstacles, until you can take the maximum height. Take your time for this.

### **Driving down off obstacles**

- Drive straight up to the edge of the sidewalk until you reach the kerb.
- Carefully Push the operating handle and allow the scooter to slowly descend from the kerb, without changing direction.

- It is not permitted to go down stairways.
- The scooter is equipped with anti-tipping wheels for extra stability and safety. Occasionally, these wheels may touch an obstacle when driving down one.
- Avoid leaving one of the back wheels on the obstacle, while the other is still on the lower level. This will cause the scooter to tip over.

### 4.3.5 Overload

- In case of overload, caused by wrong or long-term use in hot conditions, the automatic fuse will be activated, causing the scooter to stop. See 'Automatic fuse'.
- Abuse may also cause unnecessary malfunctions and defects to the scooter.


### 4.3.6 Switching off

After each ride, the scooter should be switched off completely: This prevents the batteries from unnecessarily discharging and having to be recharged too soon.

### 4.3.7 Parking

After parking, always leave the scooter switched off and 'locked'. This way, you will avoid unwanted usage and/of theft.

 After parking, pull out the contact plug from the contact, so that no one can use your

Winner without permission.

- Lock the hand brake on the front wheel as follows:
- Squeeze handle (A) completely.
- Press catch (B).
- Keep the catch pressed while releasing the handle.

After the scooter has been switched off (if the contact plug has been pulled out of the contact), the automatic parking brake will still operate, even if the batteries are removed. If the neutral position has been activated, the automatic parking brake will NOT work!

- If the scooter is stopped on a slope, the parking brake should be activated.
- Never put the scooter in neutral position if it is on a slope.

### 4.3.8 Charging the batteries

For charging, see the following documentation:

- · Regulations concerning the batteries.
- User manual for the battery charger, or the instructions on the battery charger.

The scooter has 'dry' gel batteries. These dry batteries (dry-fit) are entirely sealed and maintenance free.

With normal use, the batteries should be recharged every night.

Charging should be done as follows:

- Switch off the scooter.
- Put the plug of the charging cable into the charging connection. See 'Charging connection'.
- Put the plug of the battery charger into the wall socket of the mains power supply.
- Switch on the battery charger.

### Safety:

The operation of the scooter is switched off as soon as the plug of the charging cable is connected to the charging connection.

- The battery charger has been designed in such a way that batteries cannot be overcharged.
- When the batteries are no longer charging, the charging current will automatically drop until the charging is stopped completely: the batteries are fully charged.
- The minimum charging time for dry batteries is approximately 12 hours. Most batteries will have reached 80% of their capacity after an 8-hours' charging time.

When the batteries are charged:

- Take the plug from the charging connection.
- Switch off the battery charger.

• Take the plug of the battery charger from the wall socket of the mains power supply. The scooter is now ready for use.

Always take away the charging cable when the batteries are charged. This way, you prevent the batteries from slowly depleting themselves.

### 4.3.9 Battery charger

A good-quality battery charger is necessary to keep the batteries in optimal condition. The battery charger should meet the following requirements:

- Suitable for the charging of 2 x 12V batteries.
- The charging current should correspond with the capacity (in Ah) of the batteries. The batteries must have reached at least 80% of their capacity after an 8-hours' charging time.
- The battery charger should be doubly insulated (the charging cable does not have any mass).
- The battery charger should be suitable for automatic use: once the batteries are charged, the battery charger should switch itself off.

### 4.4 Transferring in and out

See Fig. 000000100.

Before getting in or out, or making a transfer, the scooter should be switched off electrically and the neutral position should be switched off.

Before getting in or out, or making a transfer, the armrest can be folded up, so that the seat can be turned.

# A 000000100

### Transferring in:

- Switch off the scooter.
- Move the fixation handle (A) up and turn the seat 90 degrees to the left or to the right.
   When you release the fixation handle, the seat will automatically lock: the handle will click into place.
- Now seat yourself on the seat.
- Turn back the seat as described above.
   When the seat has been turned to its original position, it will click back into place. This way, it cannot inadvertently turn while driving.

The seat is provided with armrests that can be folded up, so that you can get in and out sideways, after the armrest has been folded up.

After getting in or out, always check whether the armrests are in the proper position.

### Transferring out should be done in reverse order.

### 4.5 Pushing the scooter

In the case of calamities, a malfunction, of if the batteries have insufficient capacity to provide energy for the drive system, the scooter may also be pushed.

Pushing the scooter may also be useful if this is required because of the narrow space in which the scooter has to be parked, or if it simplifies an action:

- Switch off the scooter using the contact plug.
- Switch on the neutral. See 'Neutral position of the scooter.

If the scooter is pushed too fast, the speed will automatically slow down, or the scooter will automatically brake.

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000000110


# 4.5.1 Neutral position of the scooter

See Fig. 000000110.

To be able to push the Winner, the parking brake has to be switched off as follows:

 Move the neutral handle (A), located on the right-hand side of the hood, underneath the seat, in the RELEASE BRAKE (C) position.

This way, the automatic parking brake is switched off.

The scooter has been provided with a unique safety precaution, which prevents a speed that is too high and dangerous to be achieved while in the neutral position, i.e. with the brake function switched off. With the electronics switched on, while in neutral position, the motor brake will automatically slow down the scooter to walking speed. Next to this handle, the sticker showing the 'Neutral' indicator is mounted.

- Make sure that the neutral is in the APPLY BRAKE position when you seat yourself on the scooter.
- NEVER operate the neutral while driving.
- Do not touch the operating handle for forward or reverse driving when you operate the neutral.

The neutral handle should only be used if the scooter is to be pushed. By putting the scooter in the neutral position, the motor is disconnected mechanically, causing the parking brake of the motor not to function. That's why, after pushing, this handle should immediately be returned to the 'drive' position, which will activate the parking brake.

32

### Attention:

NEVER switch off the scooter's automatic parking brake when you are parked on a slope: the scooter will start moving down the slope as a result of gravity.

If the Winner is in the neutral position:

- It is not possible to drive electrically.
- The electronics can be switched on, but the motor cannot drive the scooter. In this situation, the electronics should be switched off.

To be able to drive the Winner electrically again, the following actions should be carried out:

- Switch off the neutral to the motor.
- Switch the scooter on.

Now, the Winner is able to drive again.

### 4.6 Reducing the Winner

The Winner can be reduced as follows:

- By removing the seat. See 'Removing the seat'.
- By folding down the steering column. See 'Folding down the steering column'.

By replacing the seat and by putting the steering column back to its operating position, the Winner is ready for immediate use again. See the descriptions concerned.

### 4.6.1 Removing the seat

See Fig. 000000120.

The seat should be removed as follows:


- Unlock the turning mechanism of the seat, using handle (A).
- Lift the seat from the socket in which it turns.
   The seat is easily lifted out if it is turned slightly in the process.


# 4.6.2 Folding down the steering column

See Fig. 000000130.

The steering column can be folded down as follows:

 Pull handle (A) down, while, at the same time, pulling down the steering column.

Use caution when folding down the steering column, to avoid getting your fingers jammed.

### 4.7 Transport

To be able to easily transport the Winner scooter in a vehicle suitable for this purpose, the Winner can be reduced. This way, it will take up a lot less space. See 'Reducing'.

The scooter can now be placed in a car. The best way to do this is by using steel planking, across which to push the scooter into the car.

The scooter may also be lifted, but this requires at least two strong persons.

- Never lift the scooter by its plastic cover.
- Considering the nature of the scooter, the user is supposed to make a transfer to a regular car seat. You may not be transported in a car or taxi while sitting on the scooter, even if the vehicle concerned has been adapted for the transport of scooters. The reason is that the Winner cannot offer the safety that is offered by standard car seats, no matter how well the Winner is fixed to the vehicle concerned.

- After the scooter has been placed inside the car, you must make sure that the scooter is not in the neutral position.
- The scooter must be secured with lashing straps to the front and the back.
- Note that the disassembled seat can cause damage as a result of car movements.
- The seat leg may be greasy.

### 4.8 Seat belt (option)

The seat belt may be mounted by the dealer to the bottom of the seat, as a standard option.

The seat belt should be mounted as follows:

- On both sides of the seat, look for the threaded mounting holes of the hinge.
- Clear these holes.
- On both sides of the hinges, mount the mounting plates of the seat belt. Fasten these
  using a bolt.

### 5 Adjustments options

### 5.1 Seat adjustments

The Winner seat offers a number of adjustment options to increase the sitting and driving comfort.

The following may be adjusted:

- A Seat height
- B Seating depth
- C Backrest
- D Lower-backrest
- E Armrests
- F Head rest
- G Steering column

### 5.1.1 Adjusting the seat height

The whole seat may be adjusted up to 10 cm in height.

This adjustment should be made by the dealer.

### 5.1.2 Adjusting the seat depth

See Fig. 000000180.

Adjust the seat depth as follows:

- Pull handle (A) up and hold.
- Slide the seat forward or back.
- Release the handle when the seat has reached the desired position. Slide the seat back or forward a little to engage the locking mechanism.

This will set the seat depth.

### 5.1.3 Backrest adjustment


See Fig. 000000140.


The backrest is continuously adjustable, from a fully horizontal to a fully vertical position. The adjustment must be done as follows:

 Pull up the handle (A) while you are on the seat

The backrest will now automatically move forward.

Press the backrest backward by leaning

back to a position that is most comfortable for you.

 After adjustment, push the handle down to secure the backrest in the position you have chosen.

### Attention:

Do not operate the handle if the seat is unoccupied, because the backrest will come forward very powerfully.

The backrest must be adjusted in such a way that your back is properly supported while driving.

If the backrest is tilted too far backwards, this will negatively influence the seating stability, especially on slopes.

### 5.1.4 Lower-backrest adjustment

See Fig. 000000150.

The lower-backrest can be adjusted to various positions, varying from more or less support in the lower-back area.

The adjustment must be done as follows:

- Seat yourself on the seat.
- Turn knob (A) in such a way that the lowerbackrest reaches the most optimal position for you.

See Fig. 000000160.


Both armrests can be folded up completely. The height and the angle of the armrest can be adjusted.


The adjustment must be done as follows:

 Turn wheel (A) to adjust the height and the angle of the armrest.

Attention: if the position of the backrest is adjusted, the position of the armrests will change along with it.

### 5.1.6 Adjusting the headrest (option)

The headrest can be adjusted as follows:

- Three height adjustments
- Tilting forward and back

### Adjusting the height of the head rest

Adjust the height of the head rest as follows:

Pull up the headrest by hand. It will engage into one of three positions with a click.

### Tilting the headrest

The head rest can be tilted as follows:

Tilt the headrest by hand until the desired position is reached.


### 5.2 Adjusting the steering column

See Fig. 000000130.

The steering column can be adjusted as follows:

- Pull handle (A) down with one hand, while at the same time pulling the steering column towards you with the other, until the most comfortable position is reached.
- Release the handle in the chosen position.
   Then, release the steering column.

If you release the steering column while still holding up the handle, the gas-spring-loaded steering column will move to its most forward position.

Use caution when adjusting the steering column, to avoid getting your fingers jammed.

### 6 Maintenance

#### 6.1 Maintenance table

Everything that is used should be maintained. This is also true for a scooter. For trouble-free use of the scooter, it should be regularly serviced by the dealer.

Below, we have indicated what needs to be checked, how often this should be done, and by whom.

Time	Description	User
Daily	Charging of the batteries after each use	X
Weekly	Checking the tyre pressure	X
	Checking for oil leakage underneath the scooter	X
Monthly	Cleaning of the scooter	X
	Cleaning of the upholstery (if necessary)	X
Every three months	• Lubricating the turning system of the seat: to do this, the seat must be removed. See 'Removing the seat'.	Х

It is recommended to have your scooter serviced by your dealer at least once a year, or, in case of intensive use, once every six months.

In principle, we recommend you let your dealer do all maintenance work. The maintenance activities you can perform yourself are indicated in the table above.

If you discover any oil leakage underneath the scooter, immediately contact your dealer. Do not drive the scooter.

#### 6.1.1 Batteries

For maintenance see the following documentation:

- Regulations concerning the batteries.
- User manual for the battery charger, or the instructions on the battery charger.

The scooter has 'dry' gel batteries. These dry batteries (dry-fit) are entirely sealed and maintenance free.

A sticker shows the connection diagram of the batteries. This sticker is placed on the inside of the cover of the battery compartment.

- Make sure that the batteries are always well charged.
- Do not use the scooter if the batteries are almost depleted. This
  may damage the batteries and you run the risk to come to an
  unintended standstill.
- The use of 'wet' batteries is not permitted. If the batteries must be replaced, dry batteries must be used.

### 6.1.2 Replacing the batteries

If the capacity of the batteries is reduced to the extend that the scooter can only make very short trips, or even none at all, the batteries are at the end of their life span. Replace them as soon as possible.

Preferably contact your dealer: he knows exactly which batteries are most suitable for your scooter.

Replace the batteries as follows:

- Switch off the scooter (remove the contact plug from the contact).
- Remove the seat. See 'Removing the seat'.
- Remove the lock (A) from the front of the battery cover.
- Remove the red and the black plugs from the battery (A).
- Loosen the taut-bands, which hold the batteries, by squeezing in the clasps (B).
- Loosen the battery clamps (C) from the cables using a spanner with a key width of 13 mm.
- Remove both batteries (D).

New batteries are placed in the reverse order, as follows:


- Put the new batteries in place.
- Fasten the batteries using the taut-bands.
- Connect the cables according to the table below:


No.	Cable colour	Text on label	Connect to
В	Black	Battery 1	Minus terminal of battery 1 (left)
С	Red	Battery 1	Plus terminal of battery 1 (left)
D	Black	Battery 2	Minus terminal of battery 2 (right)
E	Red	Battery 2	Plus terminal of battery 2 (right)

Tighten the battery clamps.

- Note the correct connection of the cables to the batteries: these are marked with labels.
- The connections may NOT be exchanged.
- Make sure that the battery terminals and battery clamps are quite clean.
- Grease the battery clamps to prevent oxidation, using acid-free Vaseline.

Replace the cover.

- As soon as the new batteries are placed, the batteries must be charged. See 'Charging the batteries'.
- Avoid metal parts from coming into contact with the battery terminals. This may lead to short circuiting, which could have very serious consequences.
- See 'Used scooters and the environment' for correct disposal of the batteries.

### 6.1.3 Cleaning the batteries

Dry batteries are maintenance free. Still, you have to pay attention to the following matters:

- Keep the batteries clean and dry: dirt and water may cause a leak, causing the capacity of the batteries to decrease.
- Cleaning the terminals: and after cleaning, grease them with acid-free Vaseline.

Never completely deplete the batteries! This can seriously damage the batteries and considerably decrease their life span.

### 6.2 Tyres

For the proper functioning of the scooter it is very important that the tyres are kept at the correct pressure.

The tyres are equipped with a car valve. You can have the tyres inflated at your local gas service station, or you can do it yourself using a foot pump. To inflate the tyres, remove the cap from the valve.

Soft tyres yield less than optimal driving for the scooter. It also costs more energy to move the scooter forward, putting a heavier load on the batteries. Moreover, tyre wear when driving on soft tyres is unnecessarily great.

For the right tyre tension, see 'Product specifications'.

Note when filling the tyres that the pressure may never exceed the maximum specified value provided in the table 'technical data', or as indicated on the sticker on the wheel concerned. When in doubt, contact your dealer / supplier.

For inspection of the tyres, see the 'maintenance table'.

### 6.2.1 Inflating the tyres

The tyres are equipped with car valves and can be inflated with a proper pump. A foot pump may be used, or you can have the tyres inflated at your local gas service station.

You may also use an ordinary bicycle pump, in combination with the reducing nipple included in the delivery.

After having pumped up the tyre, always make sure to replace the cap back on the valve, to avoid dirt and sand entering the valve.

### 6.3 Cleaning

### Removing dry dirt

Upholstery, metal parts and frame parts can usually be cleaned easily with a dry soft cloth.

### Removing mud and/or other wet dirt

Parts that are soiled with wet dirt can best be cleaned by first wiping the dirty parts with a wet sponge and then wiping them dry with a dry soft cloth.

#### **Upholstery**

Cleaning with a damp cloth and household soap. After removing the dirt the cleaned parts should be wiped dry with a soft dry cloth.

- Never use abrasive or aggressive cleansers. They can damage the scooter.
- Also do not use organic solvents such as thinner, dry-cleaning naphtha, or white spirit.
- Be careful with water near the electronic system.
- Upholstery: do not dry clean, iron or spin dry.

### 7 Troubleshooting

If your Winner will not function, even though the batteries are fully charged, check the following points before getting in touch with your dealer.

- 1. Check whether all battery clamps are securely in place.
- 2. Check whether it is in DRIVE position (not NEUTRAL).
- 3. Check if the automatic safety fuse has jumped out. If it has, push the fuse back in again.
  - If it jumps out again, contact your dealer.
- 4. The scooter has a troubleshooting diagnosis system that can be read by the dealer with special hard and software.

### 7.1 Troubleshooting table

If your scooter will not drive, or does not drive as it should, please check the list below first, to see if you can solve the problem yourself, before approaching the dealer.

Possible cause	Action	Action to be taken by
The contact plug is not properly put into the contact.	Put the contact plug into the contact properly.	Yourself
The connectors in the motor compartment are not put into the controller properly, or they are loose.	Check the connectors and make sure they are put into the controller properly.	Yourself
The battery clamps are not fastened properly.	Check the battery clamps and tighten them.	Yourself
Neutral is on.	Switch off the neutral.	Yourself
The automatic safety fuse has jumped out.	<ul> <li>See 'Automatic fuse'</li> <li>Find the cause and push the automatic fuse back in when the problem is solved.</li> </ul>	Yourself
There is a loose connection in the wiring or a connector.	Correct the connection.	Yourself
The battery voltage is too low.	<ul> <li>Check the battery voltage. See 'Battery condition meter'.</li> <li>Charge the batteries.</li> <li>Check operation of the battery charger.</li> </ul>	Yourself
The scooter drives too slowly.	Turn the speed controller up (by turning it to the right).	Yourself
Overcharge.	Check whether the scooter was frequently used on steeps slopes.	Yourself
Temperature is too high.	Allow the scooter to cool down.	Yourself

Sometimes simply switching the scooter off and back on again will correct the problem.

If you cannot solve the problem by using the list above, contact your dealer.


#### 7.2 Automatic fuse

See Fig. 000000180.

The scooter has been equipped with an automatic fuse (A).

In normal circumstances, the automatic fuse will not jump out and the scooter will be ready for use.

When the electric circuit is overloaded or shorted, the automatic fuse will switch off the whole electric installation to prevent damage to the electric installation. Overload can be caused if the motor is overheated because it is overloaded. If the automatic fuse switches off the electric circuit of the scooter, it will jump out.

The automatic fuse is located on the hood.

By pressing the automatic fuse back in again, you can try to restart the scooter. If the automatic fuse jumps out again, you need to wait awhile (15 to 20 minutes) to allow the motor to cool off.

If you are able to restart the scooter after this, you may drive on. Make sure, however, to avoid obstacles and hills.

If the automatic fuse jumps out again after the motor was allowed to cool off, there may be a short circuit. In this case, do not attempt driving, but contact your dealer.

#### 7.3 EMC malfunction

Radio waves may negatively influence the drive of an electric scooter. The scooter itself may disrupt the operation of electromagnetic fields, such as those used in alarm systems and shop circuitry.

Sources of radio waves, such as radio and TV stations, amateur radio transmitters, elevators, transmitters, stereo radios, and mobile telephones may influence electric wheelchairs and scooters.

The warnings listed below are intended to avoid the risk of unintended free-wheeling of the electric scooter, which may result in serious injuries.

1. Do not switch on manual items for personal communication, such as radios and mobile telephones, while the electric scooter is switched on.

### **Troubleshooting**

- 2. Avoid transmitters, such as radio and TV stations, in the neighbourhood. Do not get too close.
- 3. If any unintended movements should occur, or if the brake is moved to its neutral position, switch off the electric wheelchair or scooter by pulling the contact plug from the contact.
- 4. Be aware of the fact that any accessories or components that are installed, or any alterations that are made to the scooter, may increase the influence of radio waves.

#### Remarks:

There is no easy way to test the effects of radio waves on the general immunity of the electric wheelchair or scooter.

5. Notify your dealer, or the manufacturer of the scooter, of all cases of unintended movements or switching to neutral. Indicate whether there is any source of radio waves in the neighbourhood.

#### Remarks:

The immunity level of the Winner is 20 Volts/meter.

## 8 Technical specifications

## 8.1 Product specifications

Handicare B.V.:	Vossenbeemd 104
	5705 CL Helmond
	The Netherlands

### **8.1.1** Winner

Model	Winner3W / Winner 4W /	Winner
	Alpine	
Maximum user's weight	159 kg	

Description	Winner 3W	Winner 4W	Winner
			Alpine
Total length	1410 mm	-	-
Total width	673 mm	-	-
Total weight incl. batteries	153 kg	158 kg	153 kg
Total weight excl. batteries	93 kg	98 kg	93 kg
Weight of heaviest part	54 kg	59 kg	54 kg
Static stability downward Winner	> 15°	-	-
Static stability upward Winner	> 15°	-	-
Static stability sideways Winner	> 15°	-	-
Energy consumption: theoretic maximum	50 km	50km	33 km
distance Winner			
Dynamic stability on slopes Winner	> 10°	-	-
Climbing capacity for obstacles	100 mm	-	-
Maximal speed forward	15 km/h	15 km/h	9 km/h
Min brake distance from max speed	3,550 mm	-	-
Seat angle	2°	-	-
Effective sitting depth	430 mm	-	-
Effective seating width	500 mm	-	-
Sitting height from the front	455 mm-	-	-
(adjustable in steps of 25 mm (0.8")	520 mm		
Back angle	85°-180°	-	-
Back height	550 mm	-	-
Armrest height	250 mm	-	-

Description	Winner 3W	Winner 4W	Winner Alpine
Front armrest to backrest	260 mm	-	-
Minimum radius	1400 mm	1905 mm	1400 mm
Capacity to drive over obstacles	100 mm	ı	1
Ground clearance	114 mm	-	-
Turning space	<2000 mm	-	

Test data	a
Weight test dummy	159 kg

Operating force			
Operating-handle power	< 60 N		
Resetting automatic fuse	< 60 N		
Electronic switches	< 15 N		
Parking brake	< 60 N		
Inserting charging plug	< 60 N		

Technical data			
Diameter front wheel	310 x 108 mm (12.5" x 4.25")		
Diameter rear wheels	310 x 108 mm (12.5" x 4.25")		
tyre pressure wheels	276 kPa (2.76 Bar)		

Batteries	Winner 3W	Winner 4W	Winner Alpine
Maximum dimensions batteries	324 x 175 x 213 mm	-	-
Weight scooter without batteries	54 kg	59 kg	54 kg
Battery capacity	80 Ah/20h	-	-
Maximum permissible charging voltage	13.6 V	-	-
Maximum charging current	7 A	-	-
Connector type	A DIN 72311	1	-

### 8.2 Approval

The product meets the following standards:

- 1. NEN-EN 12182 (1999) Technical tools for the handicapped General standards and testing methods October 1997.
- 2. NEN-EN 12184 (1999) Electrically driven wheelchairs, scooters and associated charging equipment Standards and testing methods April 1999.
- 3. The Winner scooter is classified according to EN12184 (1999) in class C.
- 4. ISO 7176-8 Requirements and test methods for impact, static, and fatigue strength July 1998.
- 5. ISO 7176-9 climate tests for electric wheelchairs and scooters.
- 6. ISO 7176-14 Standards and test methods for controller systems of electrically driven wheelchairs 1997.
- 7. ISO 7176-16 Requirements of resistance to fire May 1997.

### 8.3 CE-Declaration

The product is in conformity with the provisions of the Medical Devices Directive and thus has CE marking.

### 8.4 Authorised service and technical support

If you have any problems or questions concerning this product, please contact your dealer. Ask us for information about your nearest dealer.

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Handicare B.V. Vossenbeemd 104 5705 CL Helmond The Netherlands