C2000

Frontantrieb /
Front-wheel drive

Heckantrieb /
Rear-wheel drive

Otto Bock®

QUALITY FOR LIFE

CE

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Instructions for Use

Table of Contents

1	General Information.....	101	4.1	Delivery	119
1.1	Preface	101	4.2	Settings.....	121
1.2	Intended Use.....	101	4.3	Initial Operation.....	121
1.3	Field of Application.....	102	5	Transportation and Storage.....	122
1.4	Technical Service	102	5.1	Folding the Backrest down	124
2	Safety	103	5.2	Additional Information.....	125
2.1	Explanation of Symbols	103	6	Operation	126
2.2	Standards and Directives	103	6.1	Setting Options	126
2.3	General Safety Instructions	103	6.1.1	Backrest	127
2.4	Safety Requirements for Transportation, Storage and Assembly.....	104	6.1.2	Armrests	128
2.5	Safety Requirements for Operation.....	106	6.1.3	Control Panel	129
2.6	Safety Requirements for Care, Maintenance and Disposal	113	6.1.4	Footrests.....	130
2.7	Requirements for the User	114	6.2	Getting Into and Out of the Power Wheelchair ...	132
2.8	Safety Functions.....	114	6.2.1	From/to the Side.....	132
2.9	Warning Symbols and Type Plates	116	6.2.2	From/to the Front	133
3	Product Description	119	6.3	Control Unit.....	133
4	Delivery and Preparation for Use.....	119	6.3.1	Control Panel	133
			6.3.2	Switching On and Off.....	135
			6.3.3	Drive Function	136
			6.3.4	“Battery Capacity” Indicator.....	138

6.3.5 Additional Electrical Functions.....	139	7.3.1 Separate LCD Monitor with Infrared.....	159
6.3.6 Drive-away Lock.....	140	7.3.2 External Mileage Counter	161
6.4 Lighting	142	7.3.3 Push-Button Module	162
6.5 Releasing and Locking the Brake	143	7.3.4 Attendant Control.....	164
6.6 Batteries	145	7.4 Other Options	164
6.6.1 Charging.....	145	7.4.1 Lap Belt	164
6.6.2 Battery Charger	146	7.4.2 Four-way Chest Strap/Chest Belt (only for Recaro® seats)	166
6.7 Seat Functions	148	7.4.3 Adapter for Headrest Mounting Kit (Standard Seat, Contour Seat)	166
6.7.1 Electric Seat Height Adjustment (Optional).....	149	7.4.4 Other Optional Add-on Components	167
6.7.2 Electric Seat Tilt (Optional).....	152	8 Malfunctions / Troubleshooting	167
6.7.3 Electric Back Angle Adjustment (Optional)	153	8.1 Warning	168
6.7.4 Recaro® Seat (Optional)	154	8.2 Error.....	168
6.7.5 Contour Seat (Optional)	156	8.3 Defect/Failure.....	173
6.7.6 Electrically Adjustable Footrest(s) (Optional)	157	9 Maintenance, Cleaning and Disinfection	173
7 Accessories.....	157	9.1 Maintenance Intervals	173
7.1 Control Panel Holders	158	9.3 Wheel Replacement	176
7.1.1 Swing-away Control Panel Holder.....	158	9.3.1 Caster Wheel Replacement	177
7.1.2 Height-Adjustable Control Panel Holder	158	9.3.2 Drive Wheel Replacement	178
7.2 Special Controls.....	158		
7.3 Other Display and Control Elements.....	159		

9.3.3 Replacing the Casing or Inner Tube of
Pneumatic Tyres 179

9.4 Replacing Defective Lamps 179

9.5 Cleaning and Disinfection..... 182

10 Technical Data 183

11 Disposal 185

12 Information on Re-use 185

13 Liability 186

14 CE Conformity 186

List of Figures

Fig. 1 Signage on the C2000 power wheelchair.. 116

Fig. 2 Main components 120

Fig. 3 Fuse housing with inserted fuse 121

Fig. 4 Inserting the cross bolt into the holder; locked
cross bolt 122

Fig. 5 Eyebolts on the driving side 124

Fig. 6 Folding the footplates up 124

Fig. 7 Folding the backrest down 125

Fig. 8 C2000 folding size..... 125

Fig. 9 Loosening the thumb screw for removal of the
side panel with armrest..... 128

Fig. 10 Adapting the armrest to the upper arm length..
128

Fig. 11 Adapting the armrests to the forearm length
129

Fig. 12 Adapting the control panel..... 130

Fig. 13 Adapting the footrest(s) to the lower leg length
131

Fig. 14 Adapting the calf pads 132

Fig. 15/ 16 Control panel; underside of control panel .. 134

Fig. 17	LCD display with all symbols	135	Fig. 36	Adapter for headrest mounting kit.....	166
Fig. 18	Drive-away lock information window	141	Fig. 37	Changing the fuse	176
Fig. 19	Front lights	142	Fig. 38	Removing the caster wheel.....	178
Fig. 20	Rear lights.....	143	Fig. 39	Disassembling the drive wheel.....	178
Fig. 21	Releasing the brake.....	144	Fig. 40	Removing the casing	179
Fig. 22	Connecting the battery charger to the power wheelchair	148	Fig. 41	Replacing the front halogen lamp	180
Fig. 23	Creep speed – Control panel display / LCD display	151	Fig. 42	Replacing the front direction indicator lamp	181
Fig. 24	C2000 with raised seat height adjustment .	152	Fig. 43	Disassembling the rear lights.....	181
Fig. 25	C2000 with tilted seat	153	Fig. 44	Replacing the direction indicator bulb	182
Fig. 26	C2000 with adjusted back angle	153			
Fig. 27	C2000 with Recaro [®] seat.....	154			
Fig. 28	Illustration of pressure points	155			
Fig. 29	Contour seat	156			
Fig. 30	Control panel holder, swing-away.....	158			
Fig. 31	Separate LCD monitor	159			
Fig. 32	Mileage counter.....	161			
Fig. 33	Push button module.....	163			
Fig. 34	Holder with attendant control.....	164			
Fig. 35	Applying the lap belt	165			

1 General Information

1.1 Preface

These instructions for use provide the user as well as his or her attendants with all the required knowledge on the design, function, operation and maintenance of the C2000 power wheelchair from Otto Bock Mobility Solutions GmbH. The instructions contain all the information that is needed to use the power wheelchair safely, to determine the possible causes of a malfunction and to help eliminate it.

Knowledge of these instructions for use is absolutely necessary for ensuring that the power wheelchair is used safely. Therefore, the user and his or her attendants must read the instructions for use thoroughly, especially the section on "Safety", before using the power wheelchair. This ensures the functionality of the power wheelchair is fully utilised.

1.2 Intended Use

The C2000 power wheelchair is intended exclusively for individual indoor and outdoor self-transportation by persons with walking impediments or walking disabilities. The C2000 power wheelchair may only be combined with the options

mentioned in these instructions for use. Otto Bock assumes no liability for combinations with medical devices and / or accessories from other manufacturers outside of the modular system. An option for the power wheelchair is to have it controlled by an attendant with the use of an attendant control.

Any other use is considered improper use. The manufacturer is not liable for any personal injury or damage to property resulting from improper use; in such cases, the user has sole liability.

The power wheelchair may only be used by properly trained persons. Training users and attendants to use the power wheelchair is required in order to protect persons from danger and to ensure that the power wheelchair is operated safely and correctly.

The operational safety of the power wheelchair can only be ensured if it is used properly in accordance with the information contained in these instructions for use. The user is ultimately responsible for accident-free operation.

1.3 Field of Application

The versatility of the C2000 power wheelchair and its modular design make it suitable for users who have walking impediments or walking disabilities due to:

- Paralysis
- Loss of limbs
- Limb defects / deformations
- Joint contractures / defects
- Other diseases

The C2000 power wheelchair was specially designed for users who are able to move independently in such a power wheelchair.

The following points should also be considered for an individual fitting:

- Body height and weight (depending on the standard version, maximum load of 160 kg or 200 kg, optionally up to 260 kg)
- Physical and mental condition
- Age of the patient
- Home and living conditions
- Environment

1.4 Technical Service

INFORMATION

Service and repairs on the power wheelchair must be carried out exclusively by experts at specialist dealers who have been authorised and trained by Otto Bock. Should any problems arise, please contact your power wheelchair supplier.

Should you have any questions or a problem that cannot be resolved despite using the instructions for use, please contact Otto Bock Customer Service (see back side of cover for address).

Otto Bock endeavours to support their customers in all respects in order to keep them satisfied with their product for a long time.

2 Safety

2.1 Explanation of Symbols

⚠ WARNING

Warnings regarding possible risks of severe accident or injury.

⚠ CAUTION

Warnings regarding possible risks of accident or injury.

NOTICE

Warnings regarding possible technical damage.

INFORMATION

Information regarding operation. Information for service personnel.

2.2 Standards and Directives

All information on safety contained in these instructions for use refers to the currently valid national laws and regulations of the European Union. In other countries, compliance with the applicable laws and national regulations is required.

In addition to the safety instructions contained in these instructions for use, the user is required to observe and comply with

the BGV (Employer's Liability Insurance Association regulations), UVV (accident prevention regulations) and environmental protection regulations. All information contained in these instructions for use must be complied with at all times without restrictions.

The power wheelchair has been constructed in accordance with state-of-the-art technology and is safe to operate. The safety of the power wheelchair has been confirmed by CE certification and the declaration of conformity.

2.3 General Safety Instructions

⚠ WARNING

Risk of suffocation. Packaging materials must be kept out of reach of children.

⚠ CAUTION

Risk of accident and injury in case of failure to follow the safety instructions. All safety instructions contained in these instructions for use and in all other applicable documents must be observed and complied with. The instructions for use must be available to the user at all times.

⚠ CAUTION

Risk of accident and injury due to improper use. The C2000 power wheelchair must only be used in accordance with the guidelines for proper use. The C2000 power wheelchair may only be used by properly trained persons. The C2000 power wheelchair may only be used to transport one person at a time.

⚠ CAUTION

Risk of burns in the proximity of fire. The back upholstery and seat cushion are flame retardant but nevertheless flammable. Therefore utmost caution is required near any sources of open flame or sparks, especially lit cigarettes.

INFORMATION

Only use original manufacturer's options. The optional components may be mounted only as described here. Failure to comply will void the warranty.

2.4 Safety Requirements for Transportation, Storage and Assembly

Transportation and Storage

⚠ WARNING

Risk of accident and injury due to improper use for transportation in motor vehicles for the disabled. For the time being, Otto Bock have not yet released the C2000 power wheelchair for use as a seat for transportation of the user in motor vehicles for the disabled! Until the power wheelchair has been released for use as a seat for transportation in motor vehicles for the disabled, wheelchair users must always transfer to the seats installed in the motor vehicle with the corresponding vehicle restraint systems! Non-observance of these instructions could lead to injury of the wheelchair user, as well as of all other passengers in the motor vehicle, in case of an accident.

You can get information on the present state of our measures from your specialist dealer or on the Otto Bock Homepage www.ottobock.com.

NOTICE

Risk of damage due to improper transportation. Hoisting devices used for transportation must have a sufficient capacity. The power wheelchair must be secured in accordance with the regulations for the transportation device used. Only attach the tensioning straps to the corresponding eyebolts and given lashing points. During transport on lifting platforms or in lifts, the wheelchair controls must be switched off. Engage the brake. Ensure that the power wheelchair is centred on the lifting platform. Ensure that no parts or components of the power wheelchair are within the danger zone.

NOTICE

Battery damage due to deep discharge. Remove the fuse for shipping or when the power wheelchair is not being used for an extended period of time.

INFORMATION

The tyres of the power wheelchair contain chemical substances that can react with other chemical substances (such as cleaning agents, acids).

Black tyres contain soot particles that can cause discolouration and black marks. A suitable surface is required when parking the wheelchair for extended periods of time.

Assembly

⚠ WARNING

Risk of accident and injury if plug or screw connections come loose. For all setting works or for folding down the backrest for transportation, screw and/or plug connections have to be loosened. This can cause uncontrolled movements of wheelchair components. To avoid injuries at pinch and shear points, carry out such works with the support of another person. Make sure to always keep parts of the body (e.g. hands or head) out of the danger zone.

⚠ CAUTION

Risk of accidents due to unsecured screw connections. When screw connections with thread locks are loosened, they must be replaced by new ones or secured with a medium-strength thread lock substance (e.g. EuroLock A24.20). After any settings on the power wheelchair are established or changed, the mounting screws and / or nuts must be retightened securely. During tightening observe torques when specified.

INFORMATION

Prior to using the power wheelchair, all the necessary mechanical adaptations (e. g. mounting special controls) and software settings (e. g. programming the control) must be made to comply with the individual requirements and abilities of the user. The settings may only be made by trained specialists who have been authorised by Otto Bock.

INFORMATION

When installing special controls, the user is obliged to have a specialist dealer give him or her detailed instructions on how to use them.

INFORMATION

The mounting position of the device for seat height adjustment and / or the seat must not be modified.

2.5 Safety Requirements for Operation**⚠ WARNING**

Risk of accidents and injury due to defective safety functions. The user is obligated to ensure that the power wheelchair and its safety functions are in safe and proper condition before every use.

The power wheelchair may only be operated if all safety functions, e. g. the automatic brakes, are functional. Brake failure can result in serious accidents with fatal injuries.

⚠ WARNING

Risk of accidents and injury due to incorrect configuration settings. Modified parameter settings in the configuration can lead to changes in driving characteristics. In particular, changes to the speed, acceleration, braking or joystick settings can lead to unexpected and therefore uncontrollable operating performance with a risk of accidents. Always test the driving characteristics of the power wheelchair after configuration / programming is complete.

Programming must only be completed by authorised personnel. Neither Otto Bock nor the control unit manufacturer are liable for damages (especially in combination with special controls) caused by programming that was not properly / professionally adapted to the abilities of the wheelchair user.

⚠ WARNING

Risk of accidents and injury due to a lack of brake functionality when the brake is disengaged. Note the lack of brake functionality when the brake is disengaged, especially when moving the power wheelchair on slopes.

⚠ WARNING

Risk of injury due to tipping during operation. The C2000 power wheelchair has been approved for ascending or descending inclines of up to 22 % for a short time. Navigating inclines above this percentage value is not permitted. The critical obstacle height of the C2000 power wheelchair is maximally 14 cm (depending on the load and driving direction). It is not permitted to cross obstacles higher than 14 cm. When driving on inclines or slopes, it is not permissible to negotiate obstacles. Operating the power wheelchair on stairs is not permitted.

⚠ WARNING

Risk of accidents due to tipping during operation. Reduce the driving speed when driving downhill (e.g. set it to speed level 1). When driving on inclines or slopes, it is not permissible to negotiate obstacles. Avoid getting into and out of the power wheelchair on inclines and slopes.

Before driving up slopes or over obstacles, a tilted seat must be lowered to its basic position, the seat height adjustment must be in its lowest position and the back must be brought into the upright position. It is advisable to slightly tilt the seat to the rear when driving downhill.

Driving over obstacles such as steps or curbs is only allowed at reduced speed (max. 3 km/h). Always approach obstacles at a right angle and cross over them without stopping.

⚠ WARNING

Risk of tipping when driving on inappropriate ground. Driving on very smooth areas (e.g. icy surfaces) or extremely rough terrain (e.g. broken rock or rubble) is not permitted.

⚠ WARNING

Risk of tipping when using lifting platforms. When using lifting platforms, lifts, busses or trains, switch off the controls of the power wheelchair and lock the brake.

⚠ WARNING

Risk of tipping due to changes in the centre of gravity. The effects of changes in the centre of gravity on the behaviour of the power wheelchair on inclines, side hills, or when crossing obstacles must be tested with the secure support of an attendant prior to initial use.

⚠ WARNING

Risk of accidents and injury due to incorrect lifting. Attendants may only lift the power wheelchair by the frame parts, but never by the footrest(s) and/or armrests.

⚠ WARNING

Risk of pinching or crushing in the electric seat adjustment range. When operating the seat height adjustment, back angle adjustment (optional) or seat tilt (optional), there are pinch and shear points in the area between the seat frame and wheelchair frame caused by the construction. All attendants must be informed about this. To avoid injuries, always make sure that no parts of the body, e.g. hands or feet, get into the danger zone, that no interfering objects, e.g. clothing, or obstacles are within the danger zone and that no unauthorised persons are present in the danger zone.

⚠ WARNING

Risk of accident and injury due to improper use of the seat height adjustment.

- The seat height adjustment function may be used only on level ground.
- When using the seat height adjustment function, no unauthorised persons may be present in the danger zone.
- There must not be any interfering objects or obstacles in the adjustment area.
- All attendants must be informed that there are pinch points caused by the construction in the area between the seat frame and wheelchair frame.

■ Users and attendants must not reach into the danger zone.

⚠ CAUTION

Risk of pinching or crushing in the pivot range of the footplate(s). Make sure your limbs are not in the danger area when folding the footplate(s) up or down.

⚠ CAUTION

Risk of falling. Do not lean the upper part of the body beyond the seat area while the seat is raised.

⚠ CAUTION

Risk of accidents due to a restricted field of vision. Keep in mind that your field of vision is restricted when driving the wheelchair with the seat in the raised position and when actuating the seat with built-in height adjustment.

⚠ CAUTION

Risk of accident when driving without experience. Driving without experience can lead to falls and other dangerous situations. Practice using the power wheelchair on level, straightforward terrain before using it.

⚠ CAUTION

Risk of accidents due to uncontrolled driving characteristics. Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunction. In this case, please contact your authorised dealer immediately. If any faults, defects or other hazards that can lead to personal injury are detected, the power wheelchair must be taken out of service immediately.

⚠ CAUTION

Risk of accidents and injury due to getting into or out of the wheelchair improperly. Switch off the wheelchair controls before getting into and out of the power wheelchair. The armrests are not capable of bearing full body weight, and therefore must not be used for getting into or out of the wheelchair.

⚠ CAUTION

Risk of accidents and injury due to rolling. A released brake may result in uncontrolled rolling of the power wheelchair. Therefore, ensure that the brake is engaged after parking the power wheelchair.

⚠ CAUTION

Risk of accidents due to uncontrolled driving characteristics. Malfunctions caused by electromagnetic fields can occur during the operation of the power wheelchair. The following particularities must be observed during operation:

- The driving characteristics of the power wheelchair can be affected by electromagnetic fields (mobile phones or other radiating devices). Switch off all mobile devices when driving.
- The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Switch off the controls whenever you do not need them.
- The power wheelchair has been tested according to EMC regulations.

⚠ CAUTION

Risk of accidents due to bad tyres. Visually inspect the tyres for sufficient tread depth and correct pressure before each use.

Whenever possible, avoid parking the wheelchair outdoors and in direct sunlight (UV light), as it causes the tyres to age prematurely. As a result, the tread surface hardens and corner pieces break out of the tread. Otto Bock recommends to replace the tyres every 2 years regardless of wear and tear.

When the wheelchair is not used for extended periods of time or if the tyres are heated strongly (e. g. near radiators or by sunlight shining through a window), permanent deformation of the tyres will result. Therefore always make sure that the C2000 power wheelchair has sufficient distance from sources of heat, move your wheelchair from time to time or jack up the wheelchair when storing it.

⚠ CAUTION

Risk of accident due to improper clothing. Otto Bock advises users to wear light-coloured clothing or clothing with reflectors during hours of darkness.

NOTICE

Damage may be caused by improper use of the seat height adjustment function. The C2000 power wheelchair may be equipped with a seat height adjustment device. The following particularities must be observed during operation:

- The seat function actuators are not intended for continuous use, but only for short-term limited operation (10 % use, 90 % idle time). General recommendation: At maximum load, 10 seconds of operation should be followed by 90 seconds of idle time. In this regard, the electric seat functions are considered independently of the driving function.
- The seat with seat height adjustment function may only be subjected to a maximum load of 160 kg / 200 kg.
- The seat height adjustment function may not be operated in case of errors or malfunctions.
- If creep speed is not activated when the seat height adjustment function is operated, please visit a specialist dealer immediately. The C2000 power wheelchair may only be used with the seat height adjustment in its lowest position until the fault is repaired.

NOTICE

Risk of damage due to excessive heat or cold. The C2000 power wheelchair may only be operated in the temperature range from –25 °C to +50 °C. It must not be operated at temperatures outside this range.

NOTICE

Risk of damage caused by overload. Depending on the standard version, the maximum load capacity for the C2000 power wheelchair is 160 kg or 200 kg. It is possible to have the C2000 power wheelchair retrofitted for a maximum load of 260 kg.

NOTICE

Breakage due to overloading the luggage carrier. The maximum load capacity for the luggage carrier is 15 kg.

NOTICE

Interference from electromagnetic fields. The power wheelchair has been tested according to EMC regulations. The following particularities must be observed during operation:

The driving characteristics of the power wheelchair can be affected by electromagnetic fields (mobile phones or other radiating devices). Switch off all mobile devices when driving.

The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Switch off the controls whenever you do not need them.

NOTICE

Damage to tyres due to excessive tyre pressure. Make sure that the tyre pressure indicated in Section 10 "Technical Data" is not exceeded.

INFORMATION

Road traffic regulations must be observed when driving in road traffic.

INFORMATION

It is strongly recommended that the user for his/her own safety always fastens the lap belt. The lap belt helps to stabilise the person sitting in the power wheelchair.

INFORMATION

Every time you switch on the control unit, it will return to the previously selected speed level. If desired, the starting speed level can be changed in the parameter settings. These settings may only be made by trained specialists who have been authorised by Otto Bock.

INFORMATION

If the driving mode is changed while driving, the power wheelchair will accelerate or decelerate.

INFORMATION

The controls of the power wheelchair are protected according to protective system IP 54 and can thus be used in bad weather (e. g. rain). The controls are approved for indoor and outdoor operation and meet the requirements regarding climate and splashed water.

INFORMATION

If the controls do not emit an error signal as soon as the brake lever is activated, this indicates a malfunction. The settings must be checked by a specialist dealer.

INFORMATION

Manoeuvring is only allowed at reduced speeds.

2.6 Safety Requirements for Care, Maintenance and Disposal

⚠ CAUTION

Risk of accidents and injury due to improper maintenance, repairs or adjustments. Only trained staff who have been authorised by Otto Bock may do maintenance work on the power wheelchair. This also applies to all repairs and settings on the brakes. Incorrect settings can lead to brake failure.

⚠ WARNING

Risk of injury due to explosive gases. Explosive gases can develop while the batteries are charging. The following safety instructions must be followed under all circumstances:

- Ensure sufficient ventilation when charging the batteries in closed rooms.
- Smoking and open flames are not permitted.
- Sparks must be avoided.
- The air vents in the trim must not be covered.

⚠ CAUTION

Risk of injury due to uncontrolled movements. The fuse must always be removed for any maintenance work where the battery cover is open.

NOTICE

Unauthorised battery replacement. The battery may only be replaced by a specialist dealer. The characteristic curve of the battery charger set at the factory corresponds to the battery provided and must not be changed. Setting the characteristic curve incorrectly can result in permanent damage to the battery.

NOTICE

Risk of damage to electronics due to water penetration. Never use a water jet or high-pressure cleaning apparatus for cleaning the C2000 power wheelchair. Water must not come into direct contact with the electronics, motor, or battery.

NOTICE

Risk of damage due to uncontrolled movements. When jacking up, secure the power wheelchair to prevent it from slipping or tilting to the side by placing a suitable base under the drive unit support.

INFORMATION

The functionality and operating safety of the power wheelchair must be verified by an authorised specialist at least once per year.

INFORMATION

Piston rods do not require lubrication. They are maintenance-free.

INFORMATION

Defective batteries must be disposed of properly in accordance with country-specific regulations.

2.7 Requirements for the User

⚠ CAUTION**Risk of accidents and injury due to improper operation.**

The power wheelchair may only be operated by a properly trained user. For this purpose, the user and any attendants must receive instruction in the use of the power wheelchair from specialists trained and authorised by Otto Bock. The user must have read and understood all the information in the instructions for use. The wheelchair may not be operated in cases of exhaustion or under the influence of alcohol or medications.

The user must not have any mental limitations which can temporarily or permanently restrict attentiveness and judgment.

2.8 Safety Functions

INFORMATION

In dangerous situations, the power wheelchair can be turned off at any time using the on/off button. When the button is pressed, the power wheelchair brakes immediately and the electric functions will be stopped.

If any malfunctions occur such as an insufficient supply of energy to the brake, the software will recognise them and trigger the emergency stop brake or reduce the speed of the power wheelchair. At the same time, a warning signal will be emitted.

INFORMATION

After each emergency stop, the controls of the power wheelchair must be turned on again. In the event of communication problems in the bus system of the controls, the system triggers an emergency stop and thus prevents any uncontrolled functions. If the driving function is still not available after switching the controls on again, unlock the brakes to activate the push mode. In this case, contact a specialist dealer as soon as possible.

2.9 Warning Symbols and Type Plates

Fig. 1 Signage on the C2000 power wheelchair

Labels	Meaning
	<p>A Type designation B Year of manufacture C Maximal load capacity (see Section “Technical Data”) D Maximal climbing ability (see Section “Technical Data”) E Maximal speed (see Section “Technical Data”) F Permissible axle load, front / rear G Permissible total weight H Symbol for separate collection of electric and electronic devices. Components of the power wheelchair and the batteries must not be disposed of like regular domestic waste. I CE marking – product safety in conformity with EC Directives J Serial number K Read the Instructions for Use prior to using the product. Observe the safety instructions in the Instructions for Use. L Manufacturer information / address M Country of manufacture</p>
	<p>Tyre pressure for front and rear wheels (in bar)</p>

Labels	Meaning
	<p>A Electric driving operation: lock motor brake B Manual pushing mode: unlock motor brake</p>
	<p>Risk of pinching. Do not reach into the danger area.</p>

Tab. 1 Signage on the C2000

3 Product Description

The C2000 power wheelchair is suitable for indoor and outdoor use. It is compactly designed and easy to manoeuvre for use indoors. Two 12 V batteries power its high-performance drive system, which makes it possible, along with spring-mounted drive wheels, to overcome obstacles (category C of EN 12184) and to ensure safe operation.

The C2000 power wheelchair is controlled by an enAble50 wheelchair control system. It includes a control panel to enter driving commands and display the current status as well as a controller that controls the drive motors and other electrical functions based on the input data. Data transmission is realised via a bus system.

Since the enAble50 can be programmed, it is possible to adapt the control unit to the personal requirements of the user; e.g. the speed, acceleration and deceleration values can all be adapted.

The C2000 power wheelchair features a chain steering.

The special features of the C2000 power wheelchair include:

- Possibility of turning the seat around using a parts kit, in order to change between front-wheel drive and rear-wheel drive.

- Easy servicing due to good and clear accessibility of all component groups
- Customisation with options and custom fabrication using modular components.

The modular design of the C2000 power wheelchair makes it possible to add other modules and devices in addition to the main components (see Fig. 2), e.g. special controls (see Section 7.2).

4 Delivery and Preparation for Use

4.1 Delivery

⚠ WARNING

Risk of accidents and injury due to incorrect settings. All later adaptation or setting work may only be carried out by authorised specialist dealers.

INFORMATION

The options included in the scope of delivery depend on the product configuration that has been selected for the C2000 power wheelchair.

Scope of delivery includes:

- Preset C2000 power wheelchair with main components (see Fig. 2)
- Battery charger
- Instructions for Use
- Optional seat functions (see Section 6.7)

Upon delivery by the specialist dealer the C2000 power wheelchair is ready for use. All settings correspond to the indications on the order form or are adjusted directly on site by the dealer. The C2000 power wheelchair is adapted to the personal requirements of the user.

The functions of the individual components can be tested by following the instructions in Section 7. Possible malfunctions are described in Section 8.


Fig. 2 Main components

- 1 Backrest
- 2 Joystick and control unit
- 3 Armrest (side panel)
- 4 Seat cushion
- 5 Fuse holder
- 6 Charging receptacle
- 7 Footrests
- 8 Luggage carrier
- 9 Drive wheel
- 10 Caster wheel
- 11 Front lighting
- 12 Rear lighting

4.2 Settings

The following tools are required to adjust the wheelchair to the requirements of the user and for maintenance work:

- Allen wrenches, sizes 3 to 6
- Open-end or combination wrenches, sizes 8, 10, 13
- Phillips screwdriver
- Screwdriver for pan head screws
- Tyre repair accessories

4.3 Initial Operation

⚠ WARNING

Risk of suffocation. Packaging materials must be kept out of reach of children.

All components must be checked for completeness and functionality prior to initial operation. Before the C2000 power wheelchair can be switched on, the 100 A fuse must be inserted into the fuse holder (see Fig. 3). The fuse holder is located laterally on the drive unit sustainer, next to the charging receptacle.

To do so, remove the fuse from the protective cover on the control panel, open the cap of the fuse housing and insert the fuse.

Make sure that the fuse is pressed into the centre of the spring contacts and that it is not at an angle.


Fig. 3 Fuse housing with inserted fuse

- 1 Fuse housing with open cap
- 2 Fuse

Initial Operation with Electric Back Angle Adjustment

⚠ CAUTION

Risk of accidents and injury due to incorrect backrest assembly. Make sure that the locking mechanism engages with the cross bolt. Verify that the cross bolt and locking mechanism are securely in place.

Folding up the backrest:

1. Fold the backrest up.
2. Insert the cross bolt at the end of the drive unit into the holder.
3. Press the lever down to have the locking mechanism engage with the cross bolt.

Fig. 4 Inserting the cross bolt into the holder; locked cross bolt

5 Transportation and Storage

⚠ WARNING

Risk of accident and injury due to improper use for transportation in motor vehicles for the disabled. For the time being, Otto Bock have not yet released the C2000 power wheelchair for use as a seat for transportation of the user in motor vehicles for the disabled! Until the power wheelchair has been released for use as a seat for transportation in motor vehicles for the disabled, wheelchair users must always transfer to the seats installed in the motor vehicle with the corresponding vehicle restraint systems!

Non-observance of these instructions could lead to injury of the wheelchair user, as well as of all other passengers in the motor vehicle, in case of an accident.

You can get information on the present state of our measures from your specialist dealer or on the Otto Bock Homepage www.ottobock.com.

⚠ WARNING

Risk of accident and injury if plug or screw connections come loose. For all setting works or for folding down the backrest for transportation, screw and/or plug connections have to be loosened. This can cause uncontrolled movements of wheelchair components. To avoid injuries at pinch and shear points, carry out such works with the support of another person. Make sure to always keep parts of the body (e.g. hands or head) out of the danger zone.

⚠ CAUTION

Risk of injury due to insufficient restraints. Secure the power wheelchair during transport in another vehicle sufficiently with tensioning straps. Only attach the tensioning straps to the corresponding eyebolts and given lashing points.

⚠ CAUTION

Risk of pinching or crushing in the pivot range of the footplate(s). Make sure your limbs are not in the danger area when folding the footplate(s) up or down.

NOTICE

Risk of damage due to falling. The maximum net weight of the C2000 power wheelchair is approx. 190 kg. Hoisting devices and transportation vehicles used for transportation must have sufficient capacity.

The fixing eyebolts on the driving side of the C2000 power wheelchair and the rims of the caster wheels serve as attachment points for tensioning straps.

Fig. 5 Eyebolts on the driving side

Before transporting the power wheelchair, switch off the control unit and engage the brake.

You can reduce the size of the power wheelchair for transportation by folding the footplate(s) up, removing the side panels and folding down the backrest.

Fig. 6 Folding the footplates up

5.1 Folding the Backrest down

Proceed as follows to fold down the backrest:

1. Remove the side panels (see Section 6.1.2).
2. Pull the release strap to the rear.
3. With your hand, fold the backrest down onto the seat surface.

Fig. 7 Folding the backrest down

Folding down the backrest with electric or mechanical back angle adjustment:

1. Remove the side panels (see Section 6.1.2).
2. Release the cross bolt at the lower end of the drive unit by pushing the lever up (see Section 4.3 "Initial Operation", Fig. 4).
3. Remove the cross bolt from the holder.
4. With your hand, fold the backrest down onto the seat surface.

Fig. 8 C2000 folding size

5.2 Additional Information

The C2000 power wheelchair must be stored in a dry place.

For transportation and storage, an ambient temperature range from $-40\text{ }^{\circ}\text{C}$ to $+65\text{ }^{\circ}\text{C}$ must be observed.

INFORMATION

If your power wheelchair is not moved for several days, permanent colour changes may occur where the wheelchair comes into contact with the surface it is standing on. Therefore a suitable surface is required when parking the wheelchair for extended periods of time.

INFORMATION

Tyres contain chemical substances that can react with other chemical substances (such as cleaning agents, acids, etc.).

Black tyres contain soot particles. They can leave black traces on the ground. For this reason, Otto Bock recommend choosing grey tyres when using the power wheelchair mainly indoors.

INFORMATION

Direct sunlight / UV light causes the tyres to age prematurely. As a result, the tread surface hardens and corner pieces break out of the tread.

INFORMATION

Avoid parking the wheelchair outdoors whenever possible. Regardless of wear and tear, replace the tyres every 2 years.

INFORMATION

Remove the fuse for shipping or when the power wheelchair is not being used for an extended period of time.

6 Operation

6.1 Setting Options

⚠ WARNING

Risk of accidents and injury due to incorrect settings. All later adaptation or setting work may only be carried out by authorised specialist dealers.

⚠ WARNING

Risk of accident and injury if plug or screw connections come loose. For all setting works or for folding down the backrest for transportation, screw and/or plug connections have to be loosened. This can cause uncontrolled movements of wheelchair components.

To avoid injuries at pinch and shear points, carry out such works with the support of another person. Make sure to always keep parts of the body (e.g. hands or head) out of the danger zone.

⚠ CAUTION

Risk of accidents due to unsecured screw connections. When screw connections with thread locks are loosened, they must be replaced by new ones or secured with a medium-strength thread lock substance (e.g. EuroLock A24.20). After any settings on the C2000 power wheelchair are established or changed, the mounting screws and/or nuts must be retightened securely. During tightening observe torques when specified.

Various adjustments can be made to the power wheelchair.

The seat height, seat width and seat angle have been set in accordance with the indications on the customer order form and may only be changed by the specialist dealer.

The rear wheel suspension, which is included in the standard equipment, must be adjusted by an authorised dealer.

The following items can be adapted by the user:

- Back angle
- Height and position of the armrests
- Lower leg length

If need be, the side panels can be removed.

6.1.1 Backrest

The backrest can be set to four different angles by pulling the release strap at the base of the backrest (see Fig. 7).

Once the backrest is in the desired position, release the strap to engage the locks in the next locking position.

6.1.2 Armrests

Removing the Side Panel with Armrest

1. Loosen the thumb screw at the lower end of the side panel holder.
2. Pull the side panel with armrest out upwards.

Fig. 9 Loosening the thumb screw for removal of the side panel with armrest

Adapting the Armrests to the Upper Arm Length

1. Loosen the Allen head screw using a 3 mm Allen wrench.
2. Move the side panel with armrest up or down along the rail into the desired position.
3. Firmly retighten the Allen head screw.

Fig. 10 Adapting the armrest to the upper arm length

Adapting the Armrests to the Forearm Length

1. Loosen the two set screws on the underside of the armrest using a 3 mm Allen wrench. On the side where the control panel is mounted, three set screws have to be loosened (see Fig. 11).
2. Move the armrest with the control panel receiver or the replacement tube along the attachment rail forwards or backwards into the desired position.
3. Firmly retighten all set screws.

Fig. 11 Adapting the armrests to the forearm length

6.1.3 Control Panel

CAUTION

Risk of accidents and injury due to uncontrolled moving of the wheelchair. Turn off the control unit of the power wheelchair before mounting or removing the control panel or changing its position.

NOTICE

Cable damage. Positioning the cable incorrectly can lead to pinching and thus damage to the cable. The cable must not be attached too tightly or too loosely. Avoid bending or squeezing the cable.

Adapting the Control Panel to the Arm Length

The control panel is attached to a rail which is located beneath the armrest.

1. Loosen the three set screws on the underside of the armrest.
2. Move the control panel with the control panel receiver along the attachment rail forwards or backwards into the desired position.

3. Firmly retighten all set screws.

Fig. 12 Adapting the control panel

INFORMATION

If the control panel rail is too long, the protruding part can be cut off with a saw.

Changing the Control Panel's Mounting Side**INFORMATION**

Only authorised personnel may change the control panel's mounting side.

The standard configuration has the control panel mounted on the right side. If desired by the user, it can also be mounted on the other side of the power wheelchair.

6.1.4 Footrests**⚠ CAUTION****Risk of pinching!**

Make sure your fingers are not in the danger area when folding the footplate(s) up or down.

The C2000 can be equipped either with a narrow or wide single-panel footrest or with two individual footrests. They are mounted to the C2000 power wheelchair with a receiver attached to rails that are already adapted in their length to the user. For the C2000 power wheelchair version with rear-wheel drive, only the narrow footplate is available. All footrest versions can be folded up.

Adapting the Footrest(s) to the Lower Leg Length

1. Using a 13 mm open-end or combination wrench and a 6 mm Allen wrench, loosen the two screws between the rails and remove them.
2. Move the footrest receiver with the footplate(s) along the rail up or down to match another fixing hole.
3. Insert the screws into the fixing holes and firmly retighten them.

Fig. 13 Adapting the footrest(s) to the lower leg length

Adapting the Calf Pads

The calf pads belong to the standard equipment of the C2000 power wheelchair. They are mounted to the rails for the attachment of the footrests. Their position can be adapted to the requirements of the user.

1. Using a 13 mm open-end or combination wrench loosen the two screws from the mounting sheet and remove them.
2. To change the height, move the receiver with the calf pad along the rail up or down to match another hole.
3. To change the depth setting, choose another hole in the mounting sheet.
4. Insert the screws into the fixing holes and firmly retighten them.

Fig. 14 Adapting the calf pads

6.2 Getting Into and Out of the Power Wheelchair

⚠ CAUTION

Risk of injury if the C2000 power wheelchair starts rolling. Switch off the wheelchair controls before getting into and out of the power wheelchair. This will automatically engage the motor brake.

⚠ CAUTION

Risk of pinching or crushing in the pivot range of the footplate(s). Make sure your limbs are not in the danger area when folding the footplate(s) up or down.

NOTICE

Risk of breakage due to overloading. Please notice that the armrests and footplate(s) are not capable of bearing full body weight and that they therefore must not be used for getting into or out of the wheelchair.

The modular design of the C2000 power wheelchair and the ease with which you can remove the side panels make it easy to get into and out of the wheelchair from / to the side or from / to the front.

Getting into and out of the wheelchair can be done by the user individually in a way that suits him or her best.

6.2.1 From / to the Side

Bring the power wheelchair as close as possible next to where the user is sitting. Depending on the side, remove the right or left side panel. The user can then slide onto the wheelchair's seat from the side. Use of a transfer board will make this even easier.

6.2.2 From / to the Front

Folding the footplate(s) up (see Fig. 6) makes it possible for the user to get into and out of the wheelchair from / to the front.

The assistance of an attendant or a transfer lifter makes it easy for the user to get into or out of the power wheelchair. Use of a rotation plate is also possible.

6.3 Control Unit

NOTICE

Risk of impairing the driving performance of the power wheelchair. The driving characteristics of the power wheelchair can be affected by electromagnetic fields (mobile phones or other radiating devices). For this reason, please switch off all mobile devices when driving.

NOTICE

Risk of damage to other devices. The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Switch off the controls whenever you do not need them.

6.3.1 Control Panel

The power wheelchair is operated using the control panel.

The control panel consists of a keypad, LCD display and joystick. The programming receptacle and two inputs for external push-buttons are located on the underside. The power wheelchair is turned on and off using the keypad; driving commands can be entered and the status of certain functions and components can be displayed.


Fig. 15/16 Control panel; underside of control panel

- | | |
|-----------------------------|------------------------------------|
| 1 Joystick | 7 Horn |
| 2 Direction indicator, left | 8 Direction indicator, right |
| 3 Warning flasher | 9 Light |
| 4 LCD display | 10 Programming receptacle |
| 5 On / Off button | 11 Inputs for extern. push-buttons |
| 6 Mode button | |

On / Off Button

The on / off button is used for switching the power wheelchair on and off and for activating the drive-away lock.

Mode Button

Briefly pressing the mode button increases the speed level. After reaching the maximum speed level, you can change back to speed level 1 by pressing the key again (1-2-3-4-5-1-...). Pressing and holding the button (at least 2 seconds) activates the “Electrical seat function” menu (see Section 6.3.5).

Horn

The horn will sound as long as the horn button is being pressed.

Warning Flasher

Pressing the warning flasher button activates / deactivates all four direction indicators.

Light

The front and rear lights are activated and deactivated by pressing the light button.

Direction Indicators

The right or left front and rear direction indicators are activated and deactivated by pressing the “direction indicator, right” or “direction indicator, left” button respectively. The direction indicators automatically go out after 20 seconds.

LCD Display

The LCD display is the communication interface between the user and the control unit. It indicates the selected speed level, the remaining battery capacity, the status of electric options and special functions as well as warnings and errors. During the start phase, all display symbols are visible.


Fig. 17 LCD display with all symbols

- | | |
|-----------------------------|-------------------------------|
| 1 Direction indicator, left | 9 Excess temperature |
| 2 Battery capacity | 10 Direction indicator, right |
| 3 Electric backrest | 11 Drive wheel brake |

- | | |
|----------------------|-----------------|
| 4 Electric seat tilt | 12 Warning |
| 5 Control panel | 13 Creep speed |
| 6 Drive-away lock | 14 Power module |
| 7 Speed level | 15 Drive motor |
| 8 Open-end wrench | 16 Light |

6.3.2 Switching On and Off

⚠ WARNING
Risk of death due to lack of brake functionality. The brake release lever must be locked when using the power wheelchair. The automatic brakes must be operational and functional.

⚠ CAUTION
Risk of accidents due to incorrect tyre pressure. Visually inspect the tyres of the power wheelchair for sufficient tread depth and sufficient tyre pressure before each use. Incorrect tyre pressure reduces the life of the tyres and impairs the driving characteristics.

The control unit of the power wheelchair is switched on and off by pressing the on/off button (see Fig. 15, Item 5). If the control is not used for a certain time, which is adjustable, the

power wheelchair is switched off automatically. It is also possible to switch off the power wheelchair with the on / off button during driving. In this case, the wheelchair brakes immediately until it stops.

INFORMATION

The parameter settings can be changed by the specialist dealer in accordance with the user's desire in what speed level or menu he / she wants the power wheelchair to be after switching on.

6.3.3 Drive Function

⚠ CAUTION

Risk of injury due to uncontrolled operation of the power wheelchair. Switch off the controls on the power wheelchair whenever you do not need them. This will prevent you from unintentionally activating the joystick.

⚠ CAUTION

Risk of injury due to tipping during operation. Observe the following safety information when operating the power wheelchair:

- Short-term operation on inclines up to a maximum of 22 %.
- When driving downhill, reduce speed in accordance with the degree of the incline.
- Do not cross differences in height of more than 14 cm.
- Do not cross steps or curbs without braking first.

The power wheelchair has been approved for ascending or descending inclines of up to 22 %. Navigating inclines above this percentage value is not permitted.

To ensure safe downhill driving, the driving speed must be reduced in accordance with the degree of incline (e.g. set it to speed level 1).

The critical obstacle height of the power wheelchair is maximally 14 cm (depending on the load and driving direction). It is not permitted to cross obstacles higher than 14 cm. Obstacles such as steps or curbs must be crossed at reduced speeds.

If there are obstacles in the travel path, it is important to drive around them with a large margin of safety. On uneven ground, the driving behaviour of the wheelchair may get out of control. Therefore the speed must always be adapted to the ground conditions.

The battery indicator and speed level are shown in the driving menu on the LCD display.

For driving use the joystick. The further the joystick is moved from the centre position, the faster the power wheelchair will travel in that direction.

The maximum speed with full deflection of the joystick depends on the selected speed level. Releasing the joystick automatically activates the brake function, which brings the power wheelchair to a halt. When standing still, the mechanical brakes are automatically activated so the power wheelchair cannot roll away.

The power wheelchair has a programmable number of speed levels (delivery condition = 5 speed levels). The "M" button (see Fig. 15, item 6) is used to increase the speed level. After reaching the maximum speed level, the control unit switches back to speed level 1.

Adapting the Driving Characteristics

WARNING

Risk of accidents and injury due to incorrect configuration settings. Modified parameter settings in the configuration can lead to changes in driving characteristics. In particular, changes to the speed, acceleration, braking or joystick settings can lead to unexpected and therefore uncontrollable operating performance with a risk of accidents. Always test the driving characteristics of the power wheelchair after configuration / programming is complete. Programming must only be completed by authorised personnel. Neither Otto Bock nor the control unit manufacturer are liable for damages (especially in combination with special controls) caused by programming that was not properly / professionally adapted to the abilities of the wheelchair user.

Only specialist dealers can use a hand programming device, which is to be connected to the control panel, to adapt the speed, acceleration and deceleration values to the individual wishes of the user.

6.3.4 “Battery Capacity” Indicator

INFORMATION

Immediately after switching the power wheelchair on, the battery indicator shows the battery capacity saved before the wheelchair was switched off the last time. The exact battery status is displayed after approximately 2 minutes.

The battery indicator on the LCD display is divided into 7 segments that show the current battery capacity.

A charge of 100 % corresponds to 7 segments on the battery symbol. If one segment disappears, this means that the battery’s capacity has decreased by approximately 14 %. At a total range of approximately 35 km over even terrain, each segment corresponds to a range of approximately 5 km when the speed is kept constant.

When the last segment flashes, the battery must be charged immediately.

If all the segments have disappeared and the battery symbol is flashing, the battery is under voltage. Since any further use will lead to battery damage, a warning signal also appears.

If all the segments of the battery symbol are flashing, this means that the battery is over voltage. Since this will lead to

battery damage, a warning symbol also appears.

The charging process is indicated by the battery segments lighting up one after the other. When the battery is charging, the driving function is locked out.

Display symbol(s)	Information
	Driving menu with speed level and battery capacity
	Low battery capacity
	Charging process with drive stop
	Battery is under voltage, with warning symbol
	Battery is over voltage, with warning symbol
	Creep speed

Display symbol(s)	Information
	Attendant control

Tab. 2 Battery indicator on the control panel

6.3.5 Additional Electrical Functions

INFORMATION

For more information on additional electrical functions, please see Section 6.7.

The other electric additional functions, such as electric back angle adjustment, seat tilt or footrests, are also activated by pressing the button “M”, which must be held down for approximately 2 seconds.

To change between the different functions, shortly press the button “M” or move the joystick to the right. The function that is currently selected is shown on the LCD display. The respective function can be extended or retracted by moving the joystick ahead or back. The electric motor moves the respective function as long as the joystick is deflected, and stops at the end positions.

Electric back angle adjustment: When the joystick is deflected ahead or back, the backrest electrically moves ahead or back accordingly. The electric motor moves the backrest as long as the joystick is deflected, and stops at the end positions.

Electric seat tilt: Deflecting the joystick to the rear will electrically tilt the seat to the rear. Deflecting the joystick to the front tilts the seat ahead until it is level (factory setting; can be adapted by the specialist dealer).

Seat height adjustment device: Deflecting the joystick ahead or back electrically raises or lowers the seat.

Electric footrests: Deflecting the joystick ahead or back raises or lowers the electric footrests.

Coupled seat adjustment: Deflecting the joystick ahead or back moves the backrest and seat forward and down or back and up simultaneously.

You can switch back to driving mode by pressing the “M” button.

Display symbol(s)	Information
	Electric back angle adjustment
	Electric seat tilt
	Coupled electric backrest and seat tilt adjustment
	Seat height adjustment
	Footrests

Tab. 3 Display of electric seat functions on the control panel

6.3.6 Drive-away Lock

INFORMATION

The drive-away lock can be activated as described in the following, only if the “drive-away lock” parameter has been set to “On”. At the factory, the drive-away lock parameter is set to “Off”.

It is also possible that the factory setting has been changed to “On” if requested with the order or by the specialist dealer. Please ask your specialist dealer about the setting selected for your wheelchair.

The control unit of the C2000 power wheelchair features an electric drive-away lock to disable the driving function.

Activation via control panel

1. With the control unit turned on, press and hold the mode button for at least 5 seconds.
2. A short beep confirms activation of the drive-away lock (if the feature that returns a beep with each push of a button is activated, then the second beep confirms activation of the drive-away lock).
3. The control unit is turned off automatically.

Activation of the drive-away lock is indicated by the key symbol on the display.

Display symbol(s)	Information
	Drive-away lock

Tab. 4 Display of the drive-away lock on the control panel

Activation via separate LCD monitor

If a special control is used, the drive-away lock is activated using the separate LCD monitor. Select the menu item "Drive-away lock" in the menu "Settings". Activation is confirmed with an acoustic signal and the control unit is switched off.

Fig. 18 Drive-away lock information window

Deactivation via control panel

1. Press the On / Off button to turn on the control.
2. Move the joystick into the front end position until you hear a beep.
3. Move the joystick into the rear end position until you hear a beep.
4. Release the joystick.
5. Another beep confirms the successful deactivation of the drive-away lock.

- The speed level and battery indicator appear on the LCD display, and the control unit is ready for driving.

The drive-away lock is deactivated and the wheelchair can be driven.

INFORMATION

If the joystick sequence is not completed correctly, the drive-away lock remains engaged. Turn the control of the power wheelchair off and on again and make another attempt at deactivating the drive-away lock.

Deactivation via separate LCD monitor

After switching the power wheelchair on, an information window appears on the LCD monitor. Deactivation is similar to "Deactivation via control panel".

6.4 Lighting

The C2000 power wheelchair is equipped with a lighting installation as standard equipment. The control panel is used to turn on the warning flashers, lights and right/left direction indicators.

Front lights

One halogen light and one direction indicator are attached to the side panels on the right and left side.

Fig. 19 Front lights

Rear lights

On the back of the power wheelchair, two back-up lamps with integrated direction indicator are installed on the frame.

Fig. 20 Rear lights

6.5 Releasing and Locking the Brake

⚠ WARNING

Risk of death due to brake malfunction. Incorrect brake settings can lead to brake failure, and therefore to serious bodily injuries and even death. Repairs and adjustments to the brake may only be carried out by authorised specialist staff.

⚠ WARNING

Risk of accidents and injury due to lack of brake functionality. When the brake is unlocked (push mode), no brake function is available. When moving the power wheelchair on an incline, the person pushing must provide the appropriate brake force.

NOTICE

Risk of damage when parking the wheelchair with the brake deactivated. Releasing the brake may result in uncontrolled rolling of the power wheelchair. Therefore, ensure that the brake is engaged after parking the power wheelchair.

If the control unit fails or there is not enough battery power, it is possible to push the power wheelchair by disengaging

the brake by means of the mechanical release mechanism. The brake release is located on the opposite side of the fuse/charging receptacle on the lateral frame.

Releasing the brake

For safety reasons, you first need to pull out the brake release bolt.

INFORMATION

When in this position, the brake must still be engaged! Should the pushing function already be active at this stage, please contact your wheelchair supplier.

INFORMATION

When the joystick is activated while the brakes are unlocked, the control system emits an error signal on the control panel. If this is not the case, there is a malfunction that must be corrected immediately by a specialist dealer.

Now the brake release lever has been released and can be pressed fully forwards.

Fig. 21 Releasing the brake

1 Brake release bolt

2 Brake release lever

In this position, the control unit will recognise that the brake has been released and automatically deactivate the driving function. As soon as the joystick is deflected, a warning will appear on the LCD display.

INFORMATION

Once the brake release lever has been released, all braking systems are deactivated.

INFORMATION

The control can be turned off while the wheelchair is being pushed.

Locking the brake

To activate the brake, push the brake release lever to the rear until the brake release bolt engages. To reactivate the driving function, switch off the wheelchair controls and then switch them on again.

6.6 Batteries**NOTICE**

Risk of property damage due to improper maintenance work on the batteries. All maintenance work as well as battery replacement may only be carried out by authorised specialist staff.

The standard version of the C2000 power wheelchair is equipped with two maintenance-free 12 V gel batteries with a capacity of 79 Ah. As an option – which is not possible with the heavy-duty version –, 12 V gel batteries with a capacity of 110 Ah can be installed. The batteries are located under the seat of the power wheelchair.

6.6.1 Charging**⚠ WARNING**

Risk of injury due to explosive gases. Explosive gases can develop while the batteries are charging. The following safety instructions must be followed under all circumstances:

- Ensure sufficient ventilation when charging the batteries in closed rooms.
- Smoking and fire are not permitted.
- Sparks must be avoided.
- Do not cover the air vents in the trim.

NOTICE

Battery damage due to deep discharge. Extended operation with low battery capacity results in deep discharge of the battery, causing damage. There is a risk that the power wheelchair may stop due to zero battery capacity and bring the user into a dangerous situation.

The range of the power wheelchair is determined by the battery capacity. Many factors have an influence on the battery capacity. In addition to the temperature, age of the batteries and amount of use, the charging cycle has a pronounced effect on the capacity and therefore on the range.

For an optimal charging frequency the following applies:

- The batteries can be charged at any time, regardless of the remaining charge.
- It takes about 10 hours until a discharged battery (lowest segment in the battery symbol is flashing) is completely charged. Subsequently, the power wheelchair can remain connected to the charger without causing any problems; the charger has a programmed recharge phase which ensures that maximum battery capacity is maintained.
- For daily use, it is recommended to connect the charger overnight so that the full capacity is available every day.
- When the power wheelchair is not used for an extended period of time, the battery will become discharged over time. It is therefore necessary to charge the battery at least once a week to maintain full battery capacity. We also recommend removing the fuse when the power wheelchair is not going to be used for an extended period of time.
- Total discharge of the batteries should be avoided.
- When the batteries are charging, the controls on the power wheelchair must be switched off to allow the charging current to be fed completely into the battery.

The following must be observed when charging the batteries:

- Only the battery charger provided by Otto Bock may be used for charging. Failure to comply with this requirement will render the warranty null and void.
- The voltage settings on the battery charger must correspond to the voltage used in your country.

6.6.2 Battery Charger

⚠ WARNING

Explosion hazard due to sparks. Always switch off the battery charger and disconnect it from the power source before disconnecting the battery.

NOTICE

Unauthorised battery replacement. The battery may only be replaced by authorised specialist staff. The characteristic curve of the battery charger set at the factory corresponds to the battery provided and must not be changed. Setting the characteristic curve incorrectly can result in permanent damage to the battery.

NOTICE

Risk of damage to the battery charger / risk of damage caused by the battery charger. Please note the following when using the charger:

- Only use a battery charger supplied by Otto Bock which has been tested and approved by Otto Bock for the respective batteries (observe information on the charger). Failure to do so can result in a battery explosion and possible impairment of health due to contact with battery acid.
- The information on the nameplate of the charger must match the country-specific voltage of the respective power supply network.
- Only use the battery charger within the specified ranges of temperature and humidity.
- Place the rubber feet of the battery charger on a level surface.
- When setting the charger up close to a window, protect it against direct sunlight.
- Keep the charger from overheating. The vent openings on the back of the charger must not be covered.
- Turn off the control unit during the charging process so that the entire charging current flows to the battery.

- Avoid dust and dirt. Only clean with a dry cloth.

The battery charger is designed for the 12 V gel batteries installed in the C2000 power wheelchair. The charger can operate along two characteristic curves to charge the respective batteries. Upon delivery of the power wheelchair, the matching characteristic curve will be set at the factory. If the battery charger is to be used with another power wheelchair or if batteries with another capacity are installed, the characteristic curve has to be checked.

To charge the batteries, proceed as follows:

1. Switch off the controls of the power wheelchair.
2. Pull back the metal slide on the battery charger's plug and connect the plug to the charging receptacle on the power wheelchair.
3. Connect the battery charger to a wall socket and turn it on. Charging will start automatically. The current charge status is indicated on the LEDs of the battery charger (see Table 5).
4. When the charging process is complete, turn off the battery charger as follows:
5. First disconnect the battery charger from the mains supply.

Pull back the metal slide on the battery charger's plug and disconnect the battery charger from the power wheelchair.

6. Turn on the controls. The power wheelchair is ready for driving.

Fig. 22 Connecting the battery charger to the power wheelchair

1 Charging receptacle

The battery charger indicates the following states:

Display	Function
Yellow LED is lit	Battery is charging
Yellow LED flashes	Battery is charged to 90 %
Green LED is lit	Battery is fully charged
Red LED is lit	Incorrect polarity (unplug the battery charger from the wall socket)
Red LED flashes	Battery is defective, or charging time has been exceeded

Tab. 5 Battery capacity indicator on the battery charger

If none of the LEDs are lit, there is no power supply.

6.7 Seat Functions

⚠ WARNING

Risk of pinching or crushing in the electric seat adjustment range. When operating the seat height adjustment, back angle adjustment (optional) or seat tilt (optional), there are pinch and shear points in the area between the seat frame and wheelchair frame caused by the construction.

All attendants must be informed about this. To avoid injuries, always make sure that no parts of the body, e. g. hands or feet, get into the danger zone, that no interfering objects, e. g. clothing, or obstacles are within the danger zone and that no unauthorised persons are present in the danger zone.

6.7.1 Electric Seat Height Adjustment (Optional)

⚠ WARNING

Risk of injury if the power wheelchair tips over. Creep speed is automatically activated when the seat height adjustment device is operated. If this is not case, this represents a malfunction and the seat height adjustment function must not be used. Contact a specialist dealer as soon as possible to eliminate the problem.

⚠ CAUTION

Risk of injury if the power wheelchair tips over. The power wheelchair can tip over if the seat height adjustment device is in the raised position and the user can fall out of the wheelchair. The following safety instructions must be followed under all circumstances:

- Only travel short distances when the seat height adjustment device is in the raised position, e. g. when manoeuvring.
- Depending on the standard version, the maximum load capacity of the seat with installed seat height adjustment device is 160 kg or 200 kg.
- Use the seat height adjustment function only with the backrest set to a vertical position.
- Apply the lap belt, if possible.
- Use the seat height adjustment function only on level ground.
- Do not change the mounting position of the seat height adjustment device.
- For retrofit or conversion of the seat height adjustment device, observe the instructions provided in the service manual.

⚠ CAUTION

Risk of pinching in the seat height adjustment range. When using the seat height adjustment function, there are shearing and crushing points in the area below the seat of the power wheelchair. To avoid injuries, do not reach into the danger area.

No unauthorised persons are permitted in the danger area. There must not be any interfering objects or obstacles in the seat adjustment area. There are crushing points caused by the construction in the area between the seat frame and wheelchair frame. Never let any parts of the body (hands, feet) or objects get under the raised seat. This risk of crushing must be explained to the attendants.

⚠ CAUTION

Risk of injury if the power wheelchair tips over. Before driving up inclines or over curbs, a tilted seat must be lowered to its basic position, the seat height adjustment device must be in its lowest position and the back must be in the upright position. Since tipping resistance is reduced in such situations, additional caution is required.

Slightly tilt the seat to the rear when driving down inclines or curbs and reduce the speed (max. 3 km/h).

⚠ CAUTION

Risk of falling. Do not lean the upper part of the body beyond the seat area while the seat is raised.

⚠ CAUTION

Risk of accidents due to a restricted field of vision. Keep in mind that your field of vision is restricted when driving the wheelchair with the seat in the raised position and when actuating the seat with built-in height adjustment.

NOTICE

Risk of damage caused by overload. Depending on the standard version, the maximum load capacity is 160 kg or 200 kg; the maximum load capacity of the heavy-duty version is 260 kg.

NOTICE

Risk of damage due to improper maintenance. Inspect the seat height adjustment device for visible signs of damage at least once a month and ensure all screw connections are tight. Also maintain correct air pressure in the tyres.

NOTICE

Risk of damage during transportation. Always lower the seat height adjustment device fully for loading or transportation.

INFORMATION

Make sure a sufficient supply of power is maintained in order to guarantee proper operation of the seat height adjustment device.

INFORMATION

The seat height adjustment device is not intended for continuous use, but only for short-term limited operation (10 % use, 90 % idle time). For every 1 hour of operation, the seat height adjustment function can be used for a maximum of 6 minutes.

As an option, the C2000 power wheelchair may be equipped with a seat height adjustment device. The seat height adjustment function allows the seat to be raised by up to 40 cm using a motor. The joystick in the “additional function” mode is used to control this function (see Section 6.3.5).

The drive function remains available even when the seat is raised. As soon as the seat leaves the lowest position, the speed is reduced (creep speed) because the power wheelchair is less tip-resistant. This is indicated by a snail symbol on the control panel or the optional, separate LCD monitor.

Fig. 23 Creep speed – Control panel display / LCD display

Fig. 24 C2000 with raised seat height adjustment

6.7.2 Electric Seat Tilt (Optional)

⚠ WARNING

Risk of injury if the power wheelchair tips over. Before driving up inclines or over curbs, a tilted seat must be lowered to its basic position, the seat height adjustment device must be in its lowest position and the back must be in the upright position.

Slightly tilt the seat to the rear when driving down inclines or curbs and reduce the speed (max. 3 km/h).

INFORMATION

For safe operation, please also observe the information in Section 2.3 “General Safety Instructions” and Section 2.5 “Safety Requirements for Operation”.

The C2000 power wheelchair can be equipped with an electric seat tilt option. The electric seat tilt adjustment lets you bring the seat into a tilted position of up to 30° for pressure relief, for example. The push button module or joystick in the “additional function” mode (see Section 6.3.5) is used to control this function. The seat can be continuously tilted backwards. The C2000 power wheelchair has an integrated feature that shifts the centre of gravity and serves to improve stability against tipping.

Fig. 25 C2000 with tilted seat

6.7.3 Electric Back Angle Adjustment (Optional)

⚠ WARNING

Risk of injury if the power wheelchair tips over. Before driving up slopes or over curbs, a tilted seat must be lowered to its basic position, the seat height adjustment must be in its lowest position and the back must be in the upright position.

Slightly tilt the seat to the rear when driving down inclines or curbs and reduce the speed (max. 3 km/h).

INFORMATION

For safe operation, please also observe the information in Section 2.3 “General Safety Instructions” and Section 2.5 “Safety Requirements for Operation”.

The backrest can be equipped with an electric back angle adjustment option. The joystick in the “additional function” mode is used to control this function (see Section 6.3.5).

The backrest can be continuously inclined backwards up to 30°.

Fig. 26 C2000 with adjusted back angle

6.7.4 Recaro® Seat (Optional)

⚠ WARNING

Risk of injury if the power wheelchair tips over. Before driving up slopes or over curbs, a tilted seat must be lowered to its basic position, the seat height adjustment must be in its lowest position and the back must be in the upright position.

Slightly tilt the seat to the rear when driving down inclines or curbs and reduce the speed (max. 3 km/h).

⚠ CAUTION

Risk of tipping due to incorrect back angle adjustment. When standing still, the backrest may be inclined to the rear up to 30° only. During driving on level ground, the backrest may be inclined to the rear up to 20° only.

INFORMATION

For safe operation, please also observe the information in Section 2.3 "General Safety Instructions" and Section 2.5 "Safety Requirements for Operation".

INFORMATION

Only authorised personnel is allowed to mount a Recaro® seat to the power wheelchair.

The C2000 power wheelchair can be equipped with various Recaro® seat models.

Fig. 27 C2000 with Recaro® seat

To adjust the back angle of a Recaro® seat model, use the turning handle on the right or left side of the backrest.

A release handle is used to fold down the backrest. It is located on the side of the backrest. By pulling the release handle upward, the backrest can be folded down to the front. If the backrest is brought into its initial position, it automatically locks in place. Verify that the backrest is securely locked in place every time it has been unlocked.

The Recaro® seat can be equipped with an electric back angle adjustment. The control in the “additional function” mode is used to control this function (see Section 6.3.5). With the Recaro® LT backrest, the lateral guide in the lumbar section can be individually adjusted by turning the hand wheel. Turning the wheel to the front brings the two lateral guides closer to each other evenly. Turning the wheel to the rear adjusts the two lateral guides away from each other.

Seat bottom

The type X and W seat bottoms have an extendable padding section at the front. The lever used to adjust the padding section is located below the seat bottom. Pull the lever up to shift the padding section into the desired position. The padding section locks in place as soon as the adjustment lever is re-

leased. The adjustment lever must audibly lock in place and return to its initial position.

Headrest

The height setting of the headrest can be changed, and the inclination is adjustable.

A second person is needed to remove the headrest. The pressure points must be detected under the cover by feel and be pressed simultaneously. Now the second person can pull the headrest up and out.


Fig. 28 Illustration of pressure points

Removing the Recaro® Seat from the Chassis

1. Pull the release strap on the right and left side beneath the seat bottom to the front. In this way, you loosen the right and left bolt that serve to lock the seat to the frame.
2. To release the posterior seat adapter from the retainer bushings, slightly tilt the seat to the rear and slide back a little.
3. Remove the Recaro® seat.

Attaching the Recaro® Seat to the Chassis

⚠ CAUTION

Risk of injury if the locking bolts are not engaged. Beim When mounting the Recaro® seat, the seat adapter must be securely attached to the retainer bushings. The anterior locking bolts must have snapped into place up to the key ring.

Place the posterior edge of the Recaro® seat on the end of the seat frame.

1. Slide the seat forward until the posterior seat adapter hangs in the retainer bushings.
2. Tilt the seat to the front until the locking bolts engage with the seat frame up to the key ring.

6.7.5 Contour Seat (Optional)

The C2000 power wheelchair can be equipped with a contour seat of different sizes. This results in greater comfort and allows for more adaptation possibilities.

INFORMATION

Only authorised personnel is allowed to mount a contour seat to the power wheelchair.

Fig. 29 Contour seat

6.7.6 Electrically Adjustable Footrest(s) (Optional)

INFORMATION

The seat function actuators are not intended for continuous use, but only for short-term limited operation (10 % use, 90 % idle time).

General recommendation: At maximum load, 10 seconds of operation should be followed by 90 seconds of idle time. In this regard, the electric seat functions are considered independently of the driving function.

INFORMATION

Only authorised personnel is allowed to mount electrically adjustable footrests to the power wheelchair.

To avoid permanent pressure loads or to ensure anti-shock positioning, the power wheelchair offers the possibility of mounting electrically adjustable footrest(s). This function can be selected in the "additional function" mode of the control (see Section 6.3.5). Use the joystick for adjustment.

The space for getting into/out of the power wheelchair can be increased by folding the footrest(s) up.

7 Accessories

⚠ CAUTION

Risk of accidents due to unsecured screw connections. When screw connections with thread locks are loosened, they must be replaced by new ones or secured with a medium-strength thread lock substance (e.g. EuroLock A24.20). After any settings on the power wheelchair are established or changed, the mounting screws and/or nuts must be retightened securely. During tightening observe torques when specified.

INFORMATION

Only original options provided by the manufacturer may be used. The optional components may be mounted only as described here. Failure to comply will void the warranty.

INFORMATION

All available optional components are contained on the order form and in the wheelchair accessories catalogue.

The C2000 power wheelchair has been designed as a modular system. Certain component groups can be exchanged and other accessory components can be added to the power wheelchair.

The enAble50 control offers the possibility to use numerous additional electrical functions as well as specially adapted input devices.

The complete range of options is listed on the order form and in the wheelchair accessories catalogue.

7.1 Control Panel Holders

7.1.1 Swing-away Control Panel Holder

To allow the user to drive the C2000 power wheelchair closer to an object or under the edge of a table, the control panel can be swung to the side with a special control panel holder (see Fig. 30).

1. Apply a little bit of pressure to press the control panel holder to the side and release the pivot element.
2. Swing the control panel holder to the side.
3. When bringing the control panel holder back to its original position, the pivot element automatically locks into place.

Fig. 30 Control panel holder, swing-away

7.1.2 Height-Adjustable Control Panel Holder

A special control panel holder can be mounted optionally, allowing the control panel to be moved downwards. To make use of this option, loosen the attachment screw (see Fig. 31) and move the control panel to the desired height.

7.2 Special Controls

Various special controls such as the sip and puff control, chin control, mini joystick, etc. can be added to the C2000 power wheelchair.

INFORMATION

Detailed information is contained in the “Special Controls” instructions for use.

7.3 Other Display and Control Elements

7.3.1 Separate LCD Monitor with Infrared

The separate LCD monitor with infrared allows all types of infrared devices to be controlled using the joystick.

Standard household and HIFI devices can be saved in the control unit using a learning mode. Switches or dimmers are offered as accessories. The mouse emulator is another additional infrared module. It allows joystick signals to be used in order to control a computer mouse.

INFORMATION

More detailed information and operating instructions are contained in a separate manual.

Fig. 31 Separate LCD monitor

Display	Function
	Driving menu with speed level and battery capacity, left direction indicator activated
	Driving menu with speed level and low battery capacity

Display	Function
	Creep speed
	Seat height adjustment (optional)
	Electric back adjustment (optional)
	Electric seat tilt (optional)
	Special functions E. g. coupled electric back and seat tilt adjustment

Display	Function
	Left footrest (optional; available from 06/2009)
	Right footrest (optional; available from 06/2009)
	Electric footrests, coupled (optional; available from 06/2009)
	Drive-away lock

Tab. 6 Information display on the LCD monitor

7.3.2 External Mileage Counter

INFORMATION

When using a separate LCD monitor, an external mileage counter is not needed, because the LCD monitor has its own mileage counter.

The external mileage counter is attached to the protective bow for the control panel. The mileage counter shows:

- Speed
- Trip mileage
- Total mileage
- Digital clock

Fig. 32 Mileage counter

Whenever the display shows km/h or m/h in the upper left corner, the speedometer is activated. Pressing the yellow button lets you see the trip mileage, total mileage, digital clock and speed one after the other.

Setting the Trip Mileage Back to Zero

Three flashing arrows on the left side of the display indicate the trip mileage counter. To set the value back to zero, press the yellow button for at least two seconds.

Setting the Time

If the digital clock is activated on the display, you may set the time. Press the yellow button for at least two seconds. The display will show 12: and 24: in turn. Pressing the button again when the display shows 24: will programme the clock to show a 24-hour display. Pressing the button when the display shows 12: will result in a 12-hour display.

Next, the hours and minutes can be set. The time digits automatically advance one after the other. Press the yellow button when the desired digits are showed on the display.

Setting the Total Mileage

When the total mileage counter is shown on the display, you can switch between km/h and m/h and enter the wheel circumference. To do so, press the yellow button for at least 2

seconds and then again press the button to select between the alternating displays (km/h and m/h). The wheel circumference can be calculated on the basis of the wheel size and can be set in the same way the clock is.

Example for a wheel size of 14 inches:

Conversion from 14 x 25.4 mm = 355.6 mm
inches to mm
 $C = d \times \pi$ 355.6 mm x 3.1416 = 1117 mm
 (wheel circumference to be set)

7.3.3 Push-Button Module

CAUTION

Risk of accidents due to uncontrolled driving characteristics. The push button module is an option for the use of the enAble50 system. Improper adaptation of the push button module to a power wheelchair with components outside the Otto Bock modular system or subsequent modifications by the user are not permitted.

The push button module lets the user select additional electrical functions directly during normal operation.

Accessible functions

Depending on the design of the power wheelchair and pro-

gramming of the push button module, up to 5 electrical functions can be controlled during normal operation (Figure 33, item 1). Examples include:

- Seat tilt
- Electric back angle adjustment
- Lift seat function
- Right / left footrest; combination of both footrests
- Special functions (combinations) => symbols S1 – S5

The symbols used on the push button module correspond to the symbols used on the LCD monitor (see 7.3.1 "Separate LCD monitor with infrared).

Push button functions

The "M" button (= mode, Figure 33, item 2) is used to advance the individual functions (1-2-3-4-5-1-2-...).

The blue LED (Figure 33, item 3) indicates the selected function.

The selected function is executed with the arrow keys (Figure 33, item 4) (function up/down).

INFORMATION

The function of the arrow keys can be adapted according to the wishes of the user (up arrow = function up or function down – depending on programming). Programming may only be completed by authorised personnel.

Three freely positioned Buddy buttons may also be connected to the push button module as options (Figure 33, item 5). They assume the function of the mode and up/down arrow keys (Figure 33, item 2/4). Symbols on the back of the push button module indicate the corresponding functions.


Figure 33 Push button module

- 1 Available electrical functions
- 2 "Function selection" mode button
- 3 "Selected function" LED indicator
- 4 Up/Down arrow keys (see information box)
- 5 Stereo plug for the connection of Buddy buttons

7.3.4 Attendant Control

For control of the C2000 power wheelchair by attendants, a separate control panel can be mounted to the backrest of the power wheelchair.

Fig. 34 Holder with attendant control

7.4 Other Options

7.4.1 Lap Belt

⚠ WARNING

Risk of injury due to improper use of lap belt.

The lap belt helps to additionally stabilise the person sitting in the power wheelchair. The lap belt must never be used as part of a restraint system for transportation of the person in motor vehicles for the disabled.

INFORMATION

It is strongly recommended that the user for his/her own safety always fastens the lap belt. The lap belt helps to stabilise the person sitting in the power wheelchair.

INFORMATION

If the C2000 power wheelchair is designed for a speed of 10 km/h and more, the lap belt is included in the standard equipment.

The power wheelchair can be equipped with an optional lap belt. The lap belt serves exclusively as a safety belt when driving the power wheelchair.

Applying the Lap Belt

1. Push the two halves of the buckle together so they engage. The belt buckle must snap in place audibly.
2. Verify that the belt has locked by trying to pulling it apart. The lap belt should not be too tight on the body. Objects caught under the belt can cause painful pressure sores.

Opening the Lap Belt

1. Press the red release button to open the belt.

Changing the belt length

The belt length can be adjusted on both sides.

1. Position both buckle halves in the middle of the body.
2. After positioning the buckle halves at a right angle in relation to the belt, their position can be varied. Any excessive belt length will be held in place by the plastic slides.

Fig. 35 Applying the lap belt

7.4.2 Four-way Chest Strap / Chest Belt (only for Recaro® seats)

⚠ WARNING

Risk of injury due to improper use of the four-way chest strap. The four-way chest strap/chest belt helps to additionally stabilise the person sitting in the power wheelchair. The four-way chest strap/chest belt must never be used as part of a restraint system for transportation of the person in motor vehicles for the disabled.

If the C2000 power wheelchair is equipped with a Recaro® seat, there is the possibility of retrofitting a four-way chest strap or chest belt.

The four-way chest strap/chest belt makes it possible to restrain the patient in the C2000 power wheelchair. The lower strap should be placed between the pelvic bone and thigh, and should not be too tight. The belt buckle is positioned at the front of the body, in the centre. The shoulder straps should be positioned evenly on both sides on the thorax, and should not be too tight.

Applying the Four-way Chest Strap / Chest Belt

1. To close the four-way chest strap, connect all parts of the belt buckle. The belt buckle must snap in place audibly.

2. Verify that the belt has locked by trying to pulling it apart. The four-way chest strap / chest belt should not be too tight on the body. Objects caught under the belt can cause painful pressure sores.

Changing the belt length

To increase the length of the belt, turn the buckle by 90° and pull. To reduce the length of the belt, pull the overhanging belt end.

7.4.3 Adapter for Headrest Mounting Kit (Standard Seat, Contour Seat)

The power wheelchair can be equipped with an optional headrest. An adapter for attachment of the headrest mounting kit is available for this purpose.


Fig. 36 Adapter for headrest mounting kit

7.4.4 Other Optional Add-on Components

- Puncture resistant tyres: Solid rubber tyres
- Seating shell interface: For the adaptation of special seating shells. Control panel holders for seating shells are also available.
- Armrest accessories: Special adapter for the armrests in our wheelchair accessories catalogue
- Joystick accessories: Fork for tetraplegics, STICK S80, soft ball, golf ball, and flexible control stick shaft
- Protective bow for control panel: Metal bow for protecting the control panel against shocks and damage
- Crutch holder
- Folding rear view mirror
- Tray
- Attachable tray top
- Pocket for mobile phone

8 Malfunctions / Troubleshooting

⚠ CAUTION

Risk of accident and injury due to uncontrolled movements of the power wheelchair. Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunction. In this case, please contact your authorised dealer immediately.

INFORMATION

Should you encounter problems while troubleshooting or if you do not manage to completely eliminate a problem by following the measures described here, please contact your specialist dealer.

Malfunctions are indicated on the control panel LCD display. Table 6 describes the individual displays with the corresponding problem sources as well as the possible causes and measures to be taken.

If malfunctions cannot be fully corrected using the procedures described here, a specialist dealer has the ability to read the exact error code with a hand-held programming device and to carry out a targeted system analysis.

All problems that have ever occurred are saved in a list and can be retrieved, e. g. in case of a general overhaul of the power wheelchair. The saved data can be used to determine future service and maintenance intervals, for example.

8.1 Warning

A warning indicates a status or malfunction of one or several components of the power wheelchair. The function of components that have no errors is not restricted. For example, if the connection between the controller and seat motor malfunctions, this error will only be indicated if the user attempts to activate the motor. However, the driving function is still available.

8.2 Error


CAUTION

Risk of injury if the power wheelchair stops abruptly. In the event of communication problems in the bus system of the controls, the system triggers an emergency stop to prevent any uncontrolled functions. Depending on the source of the malfunction, it is possible to drive the power wheelchair out of a danger zone such as road traffic after switching the control unit on again.

If the driving function is still not available after switching on the controls again, the brakes must be released (see Section 6.5) to switch over to the pushing mode. After the power wheelchair has triggered an emergency stop, contact a specialist dealer as soon as possible!

An error affects one or several functions of the C2000 power wheelchair. The system is not fully operational until the error has been corrected.

Display symbol(s)	LCD monitor display	Warning / Error	Cause	Possible corrective action
		Controller temperature warning	Overheating due to excessive load	Cool down phase
		Motor temperature warning	Overheating due to excessive load	Cool down phase
		Joystick warning	Joystick not in zero position when switching on	Bring joystick to zero position before switching on
		Hand control device fault	Defective joystick	Contact specialist dealer

Display symbol(s)	LCD monitor display	Warning / Error	Cause	Possible corrective action
		Controller fault	Defective controller	Contact specialist dealer
		Communication error (alternating flashing signal)	Faulty connection between the hand control device and controller / Defective cabling, software or hardware	Check cabling / plug connections; Contact specialist dealer
		Battery under voltage	Battery deep discharge	Charge as soon as possible
		Battery over voltage	Voltage too high (after full charge and driving downhill, for example)	Continue driving slowly

Display symbol(s)	LCD monitor display	Warning / Error	Cause	Possible corrective action
		Back angle adjustment motor fault	Faulty cabling or plug connection; defective actuator	Check cabling/plug connections
		Seat tilt motor fault		
		Seat height adjustment motor fault		
		Electric footrest motor fault (available from June 2009)		
		Drive motor fault		

Display symbol(s)	LCD monitor display	Warning / Error	Cause	Possible corrective action
		Brake fault	Open brake release / defective brake	Close brake release; check brake (e. g. the Bowden cable)
		Emergency stop	Severe fault caused by controller, hand control device, and / or drive motor malfunction	Contact specialist dealer

Tab. 7 Status and error messages

8.3 Defect / Failure

A defect represents a serious failure of a system component. A failure is the most serious type of malfunction, which causes an immediate emergency stop of the system.

Defects/failures are shown by a continuous sequential light indicator with an audible signal.

After the malfunction has been corrected, the system is re-activated by restarting it.

9.1 Maintenance Intervals

The correct function of the power wheelchair should be checked every time before using it. The items listed in table 8 must be checked by the user at the indicated intervals.

9 Maintenance, Cleaning and Disinfection

INFORMATION

A spare parts catalogue from Otto Bock is available upon request for ordering spare parts. Only spare parts and replacement parts from Otto Bock may be used. Failure to comply with this requirement will render the warranty null and void.

INFORMATION

Should you encounter problems during maintenance, contact your authorised dealer. The power wheelchair is to be checked and serviced once a year by an authorised dealer.

Component	Activity	Daily	Weekly	Monthly
Armrest and side panel	Fastening screws tightened			X
	Armrest and control panel secured		Prior to every use	
	Check armrest for damage		X	
Drive wheels	Wheels must rotate freely and without axial run out			X
	Central nut on the driving shaft tightened			X
	Check if wheel mounts are seated securely			X
	Directional stability of the entire power wheelchair		X	
Tyres	Air pressure (printed on the sidewall of the tyre)			X
	Sufficient tread depth, at least 1 mm			X
	Check for damage			X
Batteries	Check liquid / acid level (not with gel batteries)			X
Lighting	Perform visual check for damages		X	
	Test the electrical function	X		
Electronics	Control system free of errors		Prior to every use	
	Charger does not indicate any error messages on the LCD display		X	
	Check plug connections			X
Brake	Activate brake lever while control unit is switched on	X		
	Check for active brake function while the brakes are locked			X

Component	Activity	Daily	Weekly	Monthly
Footrest	Check ratchet mechanism for functionality and secure seating			X
	Check footrests for damages			X
	Visual check for scratches on the piston rod and oil leakage			X
Electric seat height adjustments	Visual inspection of all moving components, especially cabling – check for damage			X
	Check whether screw connections are tight			X
Swivelling wheels / casters	Fork seated in the receiver without play			X
	Wheels must rotate freely and without axial run out			X
	Fastening nut tightened			X
Padding and belts	Proper condition of padding			X
	No wear on the seat belts			X
	Check belt buckle for functionality		X	
Seat attachment	Check if attachment screws are seated securely			X

Tab. 8 Maintenance measures and intervals

NOTICE

When the joystick is activated while the brakes are unlocked, the control system emits an error signal on the control panel. If this is not the case, there is a malfunction that must be corrected immediately by a specialist dealer.

9.2 Changing the Fuse

The 100 A fuse is located laterally on the drive unit sustainer, next to the charging receptacle, in the fuse holder provided for this purpose (see Fig. 37).

1. Open the cap of the fuse holder.
2. Pull out the fuse.
3. Insert the new fuse into the holder. Make sure that the fuse is pressed into the centre of the spring contacts and that it is not at an angle.
4. Close the cap until you feel it snap back into place.

Fig. 37 Changing the fuse

9.3 Wheel Replacement

NOTICE

Risk of damage due to uncontrolled movements. When jacking up, secure the power wheelchair to prevent it from slipping or tilting to the side by placing a suitable base under the drive unit support.

NOTICE

Damage to tyres due to excessive tyre pressure. Make sure that the tyre pressure indicated in Section 10 "Technical Data" is not exceeded.

INFORMATION

Direct sunlight (UV light) causes the tyres to age prematurely. As a result, the tread surface hardens and corner pieces break out of the tread.

INFORMATION

Avoid parking the wheelchair outdoors whenever possible. Regardless of wear and tear, replace the tyres every 2 years.

When the wheelchair is not used for extended periods of time or if the tyres are heated strongly (e. g. near radiators or by sunlight shining through a window), permanent deformation of the tyres will result. Therefore always make sure that the power wheelchair has sufficient distance from sources of heat, move your wheelchair from time to time or jack up the wheelchair when storing it.

9.3.1 Caster Wheel Replacement

Proceed as follows:

5. Turn off the controls.
6. Jack up the power wheelchair so that the wheel you want to replace can rotate freely.
7. Loosen the four screws around the wheel hub.
8. Pull the caster wheel to the front off the wheel hub.
9. Put the new caster wheel onto the wheel hub.
10. Tighten all four screws with a torque wrench.

Properly reassemble all components upon completion of the work.

INFORMATION

When installing the wheel, tighten all four screws to a torque of **25 Nm**.

Fig. 38 Removing the caster wheel

9.3.2 Drive Wheel Replacement

Proceed as follows:

1. Turn off the controls.
2. Jack up the power wheelchair so that the wheel you want to replace can rotate freely.
3. Loosen the five cap nuts and remove them along with the washers and spring washers.
4. Pull the drive wheel to the front off the threaded bolts.

5. Put the new drive wheel onto the threaded bolts.
6. Put on the five spring washers and five washers.
7. Tighten all five cap nuts with a torque wrench.

Properly reassemble all components upon completion of the work.

INFORMATION

When installing the wheel, tighten all five cap nuts to a torque of **25 Nm**.

Fig. 39 Disassembling the drive wheel

9.3.3 Replacing the Casing or Inner Tube of Pneumatic Tyres

The rims of the wheels are two-piece rims that can be separated by removing the Allen head screws.

1. Disassemble the respective wheel.
2. Completely let the air out of the tyre.
3. Loosen all five screws that connect the two parts of the rim.
4. Lever the casing off the rim.
5. Press the valve to the inside through the rim.
6. Pull out the defective inner tube.
7. Repair the inner tube with a standard tyre repair kit or replace it with a new one.

Properly reassemble all components upon completion of the work.

Fig. 40 Removing the casing

9.4 Replacing Defective Lamps

NOTICE

Risk of damage due to moisture. Do not expose the lamps to moisture. When reinserting the pane, make sure it sits correctly in the housing and that the screws are firmly tightened.

INFORMATION

Lamp housings or lamps can be ordered from a specialist dealer.

Front lights

To change the halogen lamp in the front light, proceed as follows:

1. Fold the pane forward with slight pressure on the black locking lever at the lower end of the lamp.
2. Remove the pane with the attached halogen lamp.
3. Replace the defective halogen lamp with a new one.
4. Reinsert the pane.

Fig. 41 Replacing the front halogen lamp

To change the halogen lamp in the front direction indicator, proceed as follows:

1. Open the front direction indicator by inserting a narrow flathead screwdriver into the recess on the direction indicator housing (Fig. 42, item A) and tipping the pane down.
2. Remove the pane (Fig. 42, item B).
3. Slightly rotate the halogen lamp to release it from the lock on the lamp socket and pull it out (Fig. 42, item C).
4. For mounting, slide the halogen lamp in and rotate it into the lock.

5. Insert the lamp socket into the housing and lock the pane in place.


Fig. 42 Replacing the front direction indicator lamp

Rear lights

To change the rear direction indicators and the rear light, proceed as follows:

1. Loosen the two screws with a small flathead screwdriver and remove the protective pane.
2. Remove the defective bulb from the holder and replace it with a new one.
3. Reattach the protective pane with the two screws.

Fig. 43 Disassembling the rear lights

- 1 Direction indicator bulb 2 Rear light bulb

Fig. 44 Replacing the direction indicator bulb

9.5 Cleaning and Disinfection

NOTICE

Risk of damage to electronics due to water penetration.

In order to avoid malfunctions, the electronic components, motors and batteries must not come into contact with water when cleaning the C2000 power wheelchair.

NOTICE

Risk of damage to components of the C2000 power wheelchair. Do not use any aggressive cleansing agents, solvents or hard brushes for cleaning the C2000 power wheelchair. The power wheelchair must never be cleaned using a water jet or high-pressure cleaning apparatus.

INFORMATION

Prior to disinfection, clean the seat and back upholstery as well as the control panel and armrests.

The power wheelchair must be cleaned at regular intervals, depending on the amount of use and the degree of soiling.

The control panel, battery charger, armrests and trim can be cleaned with a damp cloth and mild cleaning solution.

Use a dry brush to clean the seat and back upholstery.

Use a damp plastic brush to clean the wheels and frame.

10 Technical Data

Dimensions and weights	
Seat width	35 – 48 cm
Seat depth	38 – 50 cm
Seat height	from 51 cm
Armrest height	24 – 36 cm
Armrest length	26 cm
Lower leg length	30 – 55 cm
Back height	55 cm
Back angle	-9/1/11/21° or 0/10/20/30°
Total width	68 cm
Total height	Dependent on the seat in use
Total length	120 cm
Turning radius	150 cm for 180° turn
Tyre size	
Drive wheel:	16"
Caster wheel:	14"
Air pressure	Front: printed on the sidewall of the tyre
Empty weight	190 kg
Max. load capacity	160 kg / 200 kg (optional 260 kg)

Corrosion protection	
Corrosion protection	Powder-coated frame
Electrical installation	
Operating voltage	24 V
Gel Batteries	
Standard:	2 x 12 V, 79 Ah (5 h)
Optional:	2 x 12 V, 110 Ah (C20), 94 Ah (C5)
Control unit:	
Model	enAble50
Operating voltage	24 V DC
Max. output current per motor	130 A
Lighting:	
Front direction indicator	H21W 12 V BAY9s
Front light	HMP 08 2.4 W; 6 V, PX13.5s
Rear direction indicator	C21W 12 V; BA15s
Rear light	C5W 6 V; BA15s
Fuse	100 A
Driving data	
Speed	6 km/h / 10 km/h / 15 km/h
Climbing ability	22 % (short-term)

Maximum obstacle height	160 kg: Caster axle ahead 10 cm 160 kg: Drive axle ahead 14 cm
Approx. distance range	60 km
Operating temperature range	-25 °C to +50 °C
Transport and storage temperature range	-40 °C to +65 °C
Battery charger	
Model	Mentzer G3-424-20 Automatic battery charger with computer-controlled characteristic curve
Power requirements	230 V ±15 %
Mains frequency	50 / 60 Hz ±4 %
Protection class	1 (protective conductor)
Charging connection	24 V DC
Nominal charging current	20 A
Residual ripple	<1 %
Characteristic charging curve	IUoU, characteristic curve analogue to DIN 41773

Primary fuse	T6 G fuse link, 3 A, not accessible from outside
Secondary circuit breaker	Electronic, reversible reverse battery protection, short-circuit-proof, idling-proof, overheating protection
Protection type	IP 21
Ambient temperature range	-10 °C to +40 °C
Display	3 LEDs
Weight	1.9 kg
Dimensions (W x H x D)	290 x 150 x 95 mm

Tab. 9 Technical data

11 Disposal

⚠ CAUTION

Risk of polluting the environment with battery acid. The batteries of the power wheelchair contain a toxic acid. They must not be disposed of with regular domestic waste and the battery acid must not enter the sewage water system or ground. The battery manufacturer's instructions printed on the batteries must be observed.

INFORMATION

If a C2000 power wheelchair is to be disposed of, all components and materials of the wheelchair must be recycled or disposed of properly.

If the power wheelchair is no longer being used, it must be disposed of properly in accordance with national environmental regulations.

Please return defective batteries to your dealer when buying new ones.

12 Information on Re-use

The power wheelchair is suitable for re-use.

Similar to second-hand machines or cars, products that are

being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle.

Based on market observations and the current state of technology, the manufacturer has calculated that the C2000 power wheelchair can be used for a period of 5 years, provided that it is used properly and that the service and maintenance instructions are observed. Periods during which the product is stored at the dealer or with the third party payer are not included in this period. It should be clearly pointed out, however, that the power wheelchair is a reliable product far beyond this defined period of time, provided that it is cared for and maintained appropriately.

In cases of re-use, the corresponding product must first be thoroughly cleaned and disinfected. Subsequently, the product must be examined by an authorised specialist to check the condition, wear and possible damages.

Any worn and damaged components as well as components which do not fit or are not suitable for the new user must be replaced.

The service manual includes a service schedule for each model, detailed information and a list of the required tools.

13 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorised dealers or by the manufacturer.

14 CE Conformity

This product meets the requirements of the 93/42/EEC guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare LP
U.S.A.-Minneapolis, Minnesota 55447
Phone +1 800 328 4058
Fax +1 800 962 2549
usa.customerservice@ottobock.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
AUS-Baulkham Hills NSW 2153
Tel. +61 2 88182800
Fax +61 2 88182898
healthcare@ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co. Ltd.
Beijing 100026 · P.R.China
Tel. +86 10 85986880
Fax +86 10 85980040
china@ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Wanchai · Hong Kong
Tel No. +852 2598 9772
Fax No. +852 2598 7886
info@ottobock.com.hk

Otto Bock HealthCare India Pvt. Ltd.
IND-Mumbai 400 071
Tel. +91 22 2520 1268
Fax +91 22 2520 1267
information@indiaottobock.com

Otto Bock Japan K. K.
Minato-Ku, J-Tokyo
Tel. +81 3 5447-1511
Fax +81 3 5447-1512
ottobock@ottobock.co.jp

Otto Bock Korea HealthCare Inc.
ROK-37-897 Seoul
Tel. +82 2 577-3831
Fax +82 2 577-3828
ottobock@korea.com

Otto Bock South East Asia Co., Ltd.
T-10900 Bangkok
Tel. +66 2 930-3030
Fax +66 2 930-3311
obsea@ottobock.co.th

Other countries

Otto Bock HealthCare GmbH
D-37115 Duderstadt
Tel. +49 5527 848-1590
Fax +49 5527 848-1676
reha-export@ottobock.de

The C2000 is covered by the following design registration:

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Otto Bock Manufacturing Königsee GmbH has been certified by the German Society for the Certification of Quality Assurance Systems (DQS) in accordance with DIN EN ISO 9001 standard, reg. no. 779 (management system)

© Otto Bock · 647 G483-D_GB-03-1001