THUASNE

MultiComfort II G2 DII concept

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Instruction manual

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Manual de utilización

912-1 Okt. 05

WWW.HELP4U.EU

WWW.HELP4U.EU

Français

Table des matières

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English

Table of contents

- 1 Introduction**
- 2 Storage and transport**
- 3 Warnings**
- 4 Description**
- 5 Technical specifications**
- 6 Assembling the bed**
 - 6.1 Assembling the bed
 - 6.2 Assembling the lifting pole
 - 6.3 Assembling the side rails
 - 6.4 Switching on the motors
 - 6.5 Cable installation
 - 6.6 Disassembling the bed
- 7 Function**
 - 7.1 Brakes
 - 7.2 Adjusting the side rails
 - 7.3 Hand control unit
 - 7.4 Hand control unit with integrated blocking device
 - 7.5 Specific warnings
 - 7.6 Lower leg rest adjustment
- 8 Maintenance instructions**
 - 8.1 Control and maintenance
 - 8.2 General condition of the framework
 - 8.3 Side rails
 - 8.4 Brakes
 - 8.5 Actuators
 - 8.6 Mains cable
 - 8.7 Other cables
 - 8.8 Hand control unit
 - 8.9 Test the general electrical system (mandatory in Germany and recommended for other countries)
 - 8.10 Protection against water
 - 8.11 Lifting pole
 - 8.12 Backup battery
 - 8.13 Other maintenance and repairing conditions
- 9 Cleaning and disinfections**
- 10 Warranty**
- 11 Disposal of the various components at the end of the bed's useful life**
- 12 What to do if the bed breaks down**

1 Introduction

Dear client,
we are pleased that you have decided yourself for this high quality product from Thuisne.

If you read the following instructions carefully and follow them this product will fulfil reliable its function a long time.

Check of delivery and set of equipment

Depending on the version it is possible that not all functions are present. Thuisne reserves the right to carry out changes without notification. Therefore, it is possible that some parts of the text are not applicable or images slightly deviates.

When you receive the product you have to assemble it first.

Please, check the delivery at receipt if the packaging is undamaged. Complain visible damages immediately.

After unpacking please check the completeness of the delivery. Use therefore figure 1 also.

The following parts are set of equipment:

- 1 x bed base head side**
- 1 x bed base foot side**
- 1 x head part**
- 1 x foot part**
- 1 x set of side rails (consists of 4 wood laths, 4 sliders each with 4 pins)**
- 1 x lifting pole**
- 1 x instructions for use**
- 1 x hand control unit with integrated blocking device**

If parts are missing or damaged, your Thuisne-distributor has to be informed immediately.

Depending on the type of delivery the packaging consists of:

- **Paper carton**
- **PE-foil**
- **Polystyrene**
- **Wood**

All these materials can be recycled and disposed properly and introduced in the recycling. Do not dispose the parts of the packaging by the house normal garbage. Take in consideration the local disposal requirements.

The components of the THUASNE bed :

1. Foot part
2. Head part
3. Bed base head end
3. Bed base foot end
5. Hand control unit

optional:

6. wooden side rails MILANO (including 4 gliders made of plastics)
7. Lifting pole
8. Triangle grip

not pictured:

- Metal side rails PARIS
- Transport and storing kit
- Extension for the bed (10 cm, to a maximum length of 220 cm)


Fig. 1: overview

2 Storage and transport

Because it is modular, the MultiComfort II G2 bed from Thuisne can easily be transported and stored. Moreover, a transport kit ensures that the bed can be moved easily when it is no longer in use.

Assembly on the transport kit

The figure below will help you to follow the instructions :

1. Dismantle the bed (see chapter 6.6)
2. Assemble the head and foot of the bed on each side of the transport kit and tighten the winged bolts. The transport kit is standing on its castors now.
3. Place both sections of the bed in the vertical tubes of the transport kit.

Your bed is now ready to be put away and transported. It can easily be moved by one person because of its castors.

Ensure that the floor is flat, stable and horizontal when you transport the bed.


Fig. 2: Assembly on transport kit (optional)

Please comply with the following storage conditions :

- Temperature : ranging from - 10 °C to + 45 °C
- Humidity : from 40 % to 70 %


Fig. 3: bed on transport kit (optional)

3 Warnings

This logo indicates a potential hazard. Pay particular attention to the hints whenever you see this sign.

Use of the care bed presupposes a thorough knowledge of the instructions for use.

The following warnings must be read before installing the bed. Please, store this instructions for use close to you.

**Before and after using the bed and each time you move or adjust it, check that all the cables are positioned in such a way that they do not trail over the floor and will not be compressed.
After use, check that all the screw joints are tightened correctly.**

The bed comprises moveable sections that allow the bed base and the height of the bed to be adjusted.
**When you are adjusting the bed, ensure that no object, animal or person is at risk of being trapped.
Once again, when adjusting the bed, check that no one is lying on the bed or close to the moveable parts.**

N.B. :The motors are very powerful. Do not allow children to play with them

Do not press the raising and lowering buttons at the same time.

**Never crush the electricity supply cables.
Never open a motor yourself and do not use the bed in the event of defective electrical components.**

Ensure that there is nothing trapping the 230 V and 24 V electricity supply cables and that they are not positioned in close proximity to the moveable section of the bed (head- or leg-raising device). Similarly, check that nothing is above the cables that could cause damage. When moving or arranging the bed, do not forget to wind the supply cable around the unit placed at the head of the bed, beneath the bed base.

Do not sit or rest on the sides of the bed when it is being adjusted electrically. You or another person shall never sit down on the bed edges if one part of the bed base is inclined or raised. Wrong use can cause damage to the actuators or to the movable structure.

Avoid placing all of the body weight on the head-raising part or on the leg part when inclined. This could damage the bed.

The bed base head end should not be loaded with more than 100 kg, the bed base foot part with not more than 50 kg by the user.

Do not forget to apply the brakes **every time** when the bed is not moved. Also check that the castor brakes are applied, especially before sitting on the bed.

When using accessories and especially mattresses to help prevent bedsores, ensure that the maximum authorised load is not exceeded.

Use only original Thuisne spare parts for repair works and maintenance. The warranty is no longer valid if other brand parts or devices, others than those recommended by Thuisne, are used.

The bed must only be used with accessories supplied by Thuisne, for instance, side rails and lifting pole etc.

If you have any doubts regarding the proper functioning of the motors, remove the plugs from the sockets so as to eliminate any risk.

The bed can also be moved with patient. Before any move of the bed, please take care that the bed base is in the lowest position and the side rails are upraised.

Move the bed slowly and do not drive over obstacles higher than 2 cm. Do not use the bed to transport patients over a long distance. Use a wheel chair or other suitable aids for this purpose.

The bed must be in the lowest position when the patient gets up or lays on bed.

In order to prevent the hand-operated control unit from being activated inadvertently, place it in a safe place or attach it to the outer section of the side rail.

The lifting pole has not to be used to lift persons but to facilitate transfer from the supine to the seated position or to alter the position. The maximum permissible load is 80 kg and the maximum permissible lateral load is 35 kg.

Do not use the side rails to lift yourself up or to carry out other movements. They are not suitable for this purpose.

Check that the bed is in the lowest position and that the side rails are fixed in the upraised position when the patient is sleeping and/or not being monitored.

If the side rails are in the highest position, check once again that they are correctly fixed into position.

The size of the used mattress has to be 200 x 90 cm and the height of mattress must not exceed the maximum height authorised for the type of side rail used (5. Technical specifications).

The following checks must always be carried out before the bed is used and at least **once a day** in order to avoid any risk to the patient :

- Are the side rails free from any visible signs of damage or deformity ?
- Are the side rails fixed together correctly, even when slight pressure is exerted ?
- Are the cables, especially the spiral cable on the mains plug, free of any visible damage ?

Regularly check the following points, **at least once a week**:

- Are the distances between the tubes or plates of the side rails always correct ? (see 8.3)
- Do all the fixing systems freely function ?
- Are all the parts free of damage, tears or visible deformities ?

Always monitor the patient's physical and mental state in order to avoid placing him/her in danger, trapping himself/herself or wounding himself/herself in another way.

Do not use electro medical devices on this care bed.

Send any complaints to the Thuasne customer service department. Please inform us of any incident occurring with the bed.

Please do not hesitate to contact Thuasne should you require a list of spare parts and useful technical information.

Please give all parts of the bed to the proper disposal. Especially the electrical parts and the battery of the main drive. We take back all the electrical parts for correct disposal.

This symbol means that the bed is intended for indoor use only.

IPX4

This symbol means that the product is protected against water splashes.

This is the symbol indicating the protection class, here type B, which means that the mains unit is fitted with a primary cut-out mechanism.

=170 Kg

The maximum permissible load is 170 kg.

Class II

"Function : ED 15% for a maximum functioning period of 2 min. Max. 5 cycles per min.":

Your care bed is fitted with an OKIN system. The latter presents a service factor of 15% for a maximum of 2 minutes of continuous functioning. Thus, when you operate the bed continuously for two minutes, you must then comply with an interval of approximately 14 minutes. Moreover, do not activate the hand-operated control unit more than 5 times per minute.

We recommend that you follow this requirement in order to increase the life-span of your care bed.

Further hints for systems with hand control unit and integrated blocking device:

Before the first use of the bed remove the fastening of the key switch at the hand control unit.

Keep the key switch ready for use and safe of access by unauthorized persons.

We recommend the purchase of a second key switch, to make sure that the blocking device can be used even when one key gets lost.

Special security device:

Your bed is equipped with a so called mains isolation (NFS). Because of a electrical switch in the mains plug it is possible to disconnect the mains voltage completely from all other electrical components. With this it is ensured that the voltage is present at the electrical components only when a button of the hand control unit is pushed.

If something does not work when the bed was not used for some days, you can activate the electrical switch mechanically by pushing the green button. Afterwards the bed drives should work again.


Figure 4

4 Description

The Multicomfort II G2 care bed is specifically designed for people over 12 years of age and more than 150 cm body height suffering from physical problems or a handicap.

This short-term healthcare bed is also suitable for bed-ridden patients cared for in their own homes as well as for traditional care in a convalescence home.

It is not suitable for hospitals.

The bed has a base that is divided into four sections.

Any patient who may use the bed must not weigh more than 135 kg (see also chapter 5 :Technical specifications) and must not be taller than 1.90 m.

Always take into account the physical and mental capacities of each patient before using the bed! There is a risk that very small patients could get their heads trapped in the side rails.

5 Technical specifications

Bed base:	90 x 200 cm
Lowest position of bed base	40 cm
Highest position of bed base	80 cm
Overall dimensions	105 x 214 cm
Idle weight	Bed: 84 kg Lifting pole: 6,4 kg Side rails: 9 kg Accessories: 0,25 kg
Maximum permitted load:	<u>Bed</u> Max. weight of the user: mattress 135 kg accessories 20 kg Total 15 kg <u>Lifting pole and triangle grip: 80 kg</u>
Maximal load of the lifting pole:	80 kg
Maximum load total	170 kg
Heaviest sections of the bed:	Weight : 20 kg
Head/foot of the bed	Dimensions: 815 x 105 x 65 mm (Mantovani) 845 x 105 x 75 mm (Haydn)
Back section :	70° (see figure 5)
Foot section :	0-30° thigh section (see figure 5) 0-20° leg section : 5 notched stages (see figure 5)
Electrical category	Class II
Noise level	Less than 65 dB
Maximum permitted mattress height:	12 to 18 cm (MILANO side rails) 12 to 19 cm (PARIS side rails)
Mattress dimensions:	90 x 200 cm


Fig. 5 : Angles of the bed base

6 Assembling the bed

The bed must be assembled, used and dismantled only by well-trained, expert personnel. Look at figure 1 and follow these instructions.

Once the bed has been assembled, follow the checking instructions described in chapter 8.

6.1 Assembling the bed

There are several ways to assemble the bed depending the experience of the assembling person, the ground and the necessary amount of physical strength. The following method is the easiest way considering the physical strength.

Pay attention while assembling, that you treat the bed parts carefully. Place them carefully on the floor, to avoid damaging the drives.

- A. Place the head part of the bed against a wall and block the two castors by applying the brakes.
- B. Stick the bed base of the head end and the head part together.
- C. Ensure that the butterfly screws are fastened tight.
- D. Position the foot part of the bed against a wall and block the two castors by applying the brakes.
- E. Stick the bed base of the foot end and the foot part together.
- F. Ensure that the butterfly screws are fastened tight.
- G. Release the brakes and push the foot and the head parts of the bed together.
- H. Push both parts completely together and tighten firmly the butterfly bolts.

**Re-check that all the screws are tight !
Before using the bed take care that all transport securing (straps and fitting) are loosened or removed.**

6.2 Assembling the lifting pole

- A Insert the lower end of the lifting pole in the bushing provided for this purpose on the outside of the bed base, at the head part of the bed. (figure 6)
- B Take care that the pin slips in the notch of the bushing.
- C Slip the loop belt of the triangle handle over the upper section of the lifting pole between the two small vertical pins.
- D The length of the triangle rein can be adapted with the buckle. Choose a length that the patient can easily reach in lying position. Check that the rein is fixed safely.

Fig. 6: lifting pole

6.3 Assembling the side rails

MILANO model (wooden side rails)

Figure 7

1. On one side, insert an axis of the side rail into the sliding channel provided for this purpose at the head and foot of the bed. Make sure that the triangular side is pointing upwards. Slide the axis upwards until you feel it click into place.
2. Then place the two side bars side by side and engage the four pins of the second axis into the ends of the bars. Check that the rounded edges of the bars and the triangular tip of the axes are pointed in the same direction.
3. Then engage the pins of the axis already affixed to the bed into the open end of the bars and guide the other axis from the bottom into the remaining free sliding channel. Also engage this axis until you feel it click into place.
4. By lifting the side rail and pressing the button recessed in the head and foot of the bed, you can unlock the side rail and lower it.

Figure 8

PARIS model (metal side rails)

Comment: when the patient is very small and mobile, he/she may slip through the space between the side rail and the head or foot base part of the bed.

1. Always position the side rail on the correct side. Follow the arrow marks on the frame to position them correctly. Side rails must always be lowered at the foot of the bed.
2. Unscrew the knurled screws (A) on the fixing clips of the side rail, so that it can be attached to the lateral tubes of the bed base. Check that the fixing clips on either side are located between the frame pillars without covering them. Then firmly tighten the knurled screws. There must be no gap in the region of the rail clips.
3. When you pull button (B) outwards (figure 9), the side rail automatically releases; you can then lower it towards the foot of the bed. When you raise it again, the lock clicks automatically into place.


Figure 9

When raising/locking the side rails in the high position or lowering/locking in the low position, make sure that no parts of the patient's body become trapped. This instruction applies particularly to the space at the head of the bed.

**Check that all of the screw connections are tight !
If the side rail is in the high position, check once again that it is firmly locked in position.**

6.4 Switching on the motors

Your care bed is fitted with an OKIN system.

Use Figure 10 to connect the actuators. The connectors are located at the main actuator. This one is mounted at the head end of the bed base.

The housing has 5 sockets marked with different symbols. Take care that the motors are connected correctly. However, there is no danger if the motors are not plugged in correctly, but they cannot be function properly.


Figure 10 shows the connections of the main actuator.

The bed is delivered with connected hand control. If the hand control cable is not connected, please assemble it as shown in figure 10.

All plugs must be fitted with a round impermeable sealing ring, also known as a O-ring. If this sealing is missing, there will be inadequate protection against humidity.

Check that each plug is inserted firmly that each of these O-ring slides into the cover of the corresponding socket, thus ensuring correct plug connection.

1. hand control (M)
2. back rest, blue (I)
3. leg rest, green (II)
4. height adjustment head part, black (III)
5. height adjustment foot part, black (IV)

IMPORTANT: A wrong connection can cause serious damage. Please pay attention to the correct position of each plug as pictured below.

Fig. 11: Correct connection of drives and hand control unit

NO!

Only use the below mentioned combinations of hand control and main drive. Wrong connections may cause irreparable damage to the components:

Main drive (orange battery lid)	with	Hand control (with key switch)
1.39.000.191.30		1.42.500.001.30
1.39.000.147.30		
1.39.000.186.30		1.69.500.001.21

6.5 Cable installation

The hand control unit is factory-assembled clinked to the anti-traction system on the outside of the main control. If this is not the case, please fix the cable of the hand control unit as shown in figure 12.

Fig. 12: traction control for hand control unit cable

The two coiled cables used for height adjustment must be positioned diagonally and freely for connection to the main actuator (figure 13).

Fig. 13: cable from drive for height adjustment

The network cable is lead ex factory downwards and is coiled one time around over the main drive (figure 14) and is then led towards the traction relief at the head part.

Fig. 14: laying network cable

Make sure, that the traction relief applied on the network cable (figure 15) is lead into the fixation on the frame (attached ex factory).

Fig. 15: traction relief network cable

The network cable has to be coiled around the u-shaped holder in that way that the free end between the network plug and the bed is not lying on the floor.

If the bed is moved you have to pay attention that the cable is not coiled to much, to avoid to high traction.

Fig. 16: cable holder

From there lead the cable between the two cross tubes (figure 17) at the head section to the main connection. You minimize the risk to cross the main cables while moving the bed.

Figure 17

After finalisation of all tasks, check that, when moving through all bed adjustments the cables are not pulled too tightly or squeezed, not get rubbed or trailing over the floor.

6.6 Disassembling the bed

1. Move the bed base and head- and leg-raising devices into the lowest position, namely the horizontal position.
2. Remove the plug from the mains. Undo all electrical connections (adjustable height mechanism and the leg-raising devices). Hand control unit and back rest device may stay plugged in. Wind the mains connection cable to the U-shaped support located on the head of the bed for this particular purpose.
3. Then disassemble the bed by following the assembly instructions in reverse order.

While disassembling and transporting the bed moveable parts should be secured (e.g. with straps) to avoid falling down and being damaged.

7 Function

The assembly and use of the bed normally require the following conditions to be met :

- Wall socket 240 V (UK) or 230 V (IEC), 50 Hz ;
- Flat, stable and horizontal floor in a dry room ;
- An available surface area of at least 2.5 m x 1.5 m, with a height of 2.2 m.

7.1 Brakes

Do not forget to block the castors by engaging the brakes prior to use.

Fig. 18: castor, not blocked

Push down for fixation

Fig. 19: castor blocked with applied brake

Push forward to loose

7.2 Adjusting the side rails


Figure 20

MILANO :

To lower the side rail, lift it up slightly and press the unlocking button (A). Then fully lower the side rail.

Proceed in the same way at the other end of the bed.

PARIS:

When you pull down button (B), the side rail automatically releases and you can lower it towards the foot end of the bed. When you raise it again, the lock clicks automatically into place.


Figure 21

If the side rail is raised to the highest position, always make sure that it is safely locked.

When raising/locking the side rails in the high position or lowering/locking in the low position, make sure that no parts of the patient's body become trapped. This instruction applies particularly to the space at the head of the bed.

Do not let go the side rails while you lower them, they may move down too fast. This would cause a risk of hurting the patient and/or damage the side rails.

7.3 Hand control unit

The electrical control unit supplied with the bed presents a service factor of 15 % for not more than 2 minutes of continuous operation. This must then be followed by an interval of at least 14 minutes. Moreover, the hand control must not be activated more than 5 times per minute. We advise you to comply with this instruction in order to increase the life-span of your Thuisne care bed.

Do not operate several keys at the same time, to avoid to load the system and to cause damage in consequence.

When you adjust the head or leg part of the bed base, check that the patient's limbs (arm or leg) do not extend over the edge of the bed or the patient to avoid an injury.

Adjusting the head part

Press buttons 1 (upwards) and 2 (downwards)

Adjusting the foot part of the bed base

Press buttons 3 (upwards) and 4 (downwards)

Adjusting the height

Press buttons 5 (upwards) and 6 (downwards)

Fig. 22: hand control

7.4 Blocking device

Your bed has an integrated blocking device. With this equipment all functions of the bed can be blocked centrally.

These security equipment is placed on the back side of the hand control.

With the removable key switch you can choose the desired function.

Factory-assembled the key switch is connected to the hand control. Remove the fastening to make undesired changes more difficult for unauthorized persons.

We recommend to have a second key switch ready for use, to ensure that the blocking device can be operated in case of a lost.


Fig. 23: Integrated blocking device

You must always deactivate the bed functions if you are uncertain as to the patient's ability to operate the bed control system due to his/her physical or mental condition.

7.5 Specific warnings

By a drop out of the electric tension the bed can be brought to the lowest and horizontal position by the use of the backup battery. This works only, when the battery is charged, the integrated is deactivated and all electrical components work without any fault.

Attention: This function has not to be employed for regularly use, but especially in case of emergency. Afterwards the backup battery has to be renewed immediately (see chapter 8.12).

If the motor of the head end bed base dropped out, it can be manually placed in the horizontal position. To do this, the motor must be detached from the bed base and placed on a support. Two people are needed to carry out this manoeuvre – one person to lift the bed base and the second to pull one of the attachment axes (fold back the loops and pull upwards). The head-raising device can then be placed in the horizontal position.

If you wish to re-connect the actuator to the bed base, lower the actuator using the hand control and re-connect it to the bed base with the pins.

Fig. 24: fixing of the motors

Check in particular that the loops are once again positioned so that the axes cannot come out of their receptacle.

When you disconnect the plug from the mains, do not forget to wind this cable around the U-shaped support provided for this purpose at the head of the bed (see figure) so as not to damage the cable in any way.

7.6 Lower leg rest adjustment

The bed is equipped with an adjustable leg rest for the lower legs. This part may be placed in several positions by bringing it manually in the wanted position. By bringing the leg rest entirely electric down, the next time the fowler (sloping) position will be taken automatically. Exception on this is the horizontal position: this remains until the part for the lower legs is further pulled upwards resulting back to fowler position.

8 Maintenance instructions

8.1 Control and maintenance

To avoid any accidents and to increase the life-span of the Thuisne bed, a thorough check procedure should be conducted at least once a year or each time the bed is re-installed (consult local legislation, which may impose controls).

For this purpose, use the enclosed protocol. The check procedures described must be carried out only by suitably specialised personnel who have received corresponding training or else under their supervision.

If, when carrying out the afore-mentioned tasks, you have the slightest doubt concerning the safety or correct function of just one section of the bed, you must not, under any circumstances, use it. Contact your supplier or the manufacturer.

The following tasks must be carried out in all cases :

8.2 General condition of the framework

Good general condition of the framework : Check that the platform, head and foot, backrest, seat, foot rest, bracket and castors are not twisted.

The following tests must always be carried out :

- Do all the mobile parts move easily ?
- Are there any parts that hang loosely or rattle when the bed is moved ?
- Are there any scraping or grinding noises when the bed is moved ?
- Are there any rusted parts ?
- Are all the motor supports intact and free of any other visible signs of deterioration ?
- Can the bed base be adjusted at all angles and in all positions ? Measure it !
- Can the minimal and maximum heights of the bed base be reached ? Measure it !
- Place the bed rest in its lowest position, load it evenly with the maximum weight and place it in its highest position. Calculate the length of time needed by the motors at the head and foot of the bed to reach this position. The difference must not exceed 3 seconds.

8.3 Side rails

- Check that these have not been altered. Check that the locking system is working properly. There must not be a gap of more than 10 mm.
- The deviation between two tubes (PARIS) or two plates (MILANO) must not have changed and must be less than 120 mm. The rails must function easily without any significant friction or slackness.

8.4 Brakes

- Check that the brakes and the locking system are working properly. The wheels must turn freely when the brakes are unlocked.

8.5 Actuators

- Check that their range is complete (400 mm in the up/down direction, 110 mm for the back rest, 3° for the foot rest).
- Check the noise level (which must be less than 65 dB), speed and resistance.
- Each actuator must stop at the end of its range and must not exceed the last notch.
- When bringing high the bed the maximum secondary intensity value must be less than 1,7 A.

- The plastic covers for the actuators must not be damaged.
- Check that the transformer housing is fixed firmly to the motor on the head-raising device..

8.6 Mains cable

- Check that the main cable is not sectioned and that the protective plastic has not been broken.
- Also check the condition and correct functioning of the anti-traction system on the supply cable (see figure 15) and the cord connector.

8.7 Other cables

- Check that they are in a good condition. The plastic surrounding must not be torn or crushed.
- Check the seals on the connecting plugs.

8.8 Hand control unit

Please check that

- All keys are moving easily and do not stick
- The box is without any tear or other damage
- The cable is not kinked or damaged
- The plug or the kink protection cover are not torned.

The integrated blocking device:

The function of the security equipment have to be controlled. Please proceed in the following way:

1. put the switch in position 2 (lock symbol open, see figure 23)
2. move all inclining of the bed in the lightly upraised position
3. put the switch in position 1 (lock symbol close, see figure 23)
4. by pressing each adjustment button the bed must not move
5. bring the key switch on the back side of the hand control unit in testing position 1
6. while pressing each adjustment button the bed must not move
7. bring the key switch on the back side of the hand control unit in testing position 2
8. while pressing each adjustment button the bed must not move

8.9 Test the general electrical system (mandatory in Germany and recommended for other countries)¹

Checking of the electrical system is mandatory. This must at least include examination of the emergency fault current of the control unit and determination of insulation resistance². This check procedure must be implemented by an electrician or electro-technician using appropriate measuring and test methods.

Following conditions have to be fulfilled:

Emergency fault current $\leq 0,100 \text{ mA}$

Insulation resistance $\geq 7 \text{ MOhm}$

During this check procedure, the green button on the plug and the corresponding key on the hand-operated control unit must now be pressed.

Measurements must be recorded between :

- the control unit and the bed framework
- the control unit and the hand-operated control unit.

If the measured values exceed the values indicated above or if they deviate by more than 50% from the first / previous measurement, the electrical system must be taken out of use and the system must be inspected by a Thuisne distributor.

¹ In accordance with German law VBG4 / BGV A 2 (only for Germany)

² In accordance with DIN VDE 0751

8.10 Protection against water

Check that the cover is in a not damaged and free of cracks. Inspect to ascertain that the seals are still present and in perfect condition.

8.11 Lifting pole

Check that the bushing and the inner compartment (cover plate at the upper end of the head) are free from any corrosion and cracks. Check that the handle and the fixation belt are free from deformations, damages of the buckle or cracks.

8.12 Backup battery

The backup-battery is used for the mains isolation. In case of mains supply drop out the bed can be brought in the lowest position and the bed base can be adjusted flat with it. Therefore the bed is working without mains supply but this leads very soon to the discharge of the battery.

When the backup battery is discharged, it has to be replaced by a new battery of type 9V-6LR61 (Type Alkaline). Proceed as described below:

Open the (orange) coloured battery case at the main actuator (2 cross head screws). Be careful to avoid damaging the sealing ring. Pull the old battery out, remove the connector. For mounting the new battery follow the instruction the other way round. Check that the sealing ring is not damaged and mounted correctly.

Fig. 25: exchanging the backup battery

The backup-battery should be replaced minimum every two years.

8.13 Other maintenance and repairing conditions

Furthermore, for safety reasons, the following measurements must be recorded every five years :

- Change of the lifting pool ;
- Careful check of all actuators (no evidence of cable and connection deterioration, noise level).

9 Cleaning and disinfections

Before cleaning and disinfecting, pay attention to the following measures: Check that all cables and connections are intact, to prevent liquid to enter. Damaged parts have to be replaced.

Check that the mains plug has been disconnected before cleaning the bed. Please position it to avoid that the contacts get wet. After finishing the cleaning, the plug has to be dried completely.

Make sure that all actuators and the hand control unit are connected, to prevent humidity entering the plugs and sockets to avoid defective function or even short circuits. On every plug have

sealing rings present and intact. All sockets have to be closed with the matching plug or a covering cap.

**Never use a high-pressure cleaner to clean this bed.
Do not splash the bed with water.**

Although the electrical elements are protected against splashing, clean them only with a damp cloth and if necessary with a mild non-abrasive detergent.

Before using the bed again all parts have to be perfectly dry again.

You can clean the metal parts of the bed by hand with water and some soap (for example mild detergent or car shampoo)

The wooden parts wipe with a damp cloth and wipe them dry immediately. Use suitable wood care-detergent for maintenance.

Before re-using the bed, spray the bed moderate with a mild commercial disinfectant. For plastic and wooden parts do not use high percentage alcohols, ether, ketenes, ester, high concentrated acids and aromatic or chlorinated hydrocarbons.

The wooden parts have to be treated with a disinfectant with a fast effect. Sprinkle the wood and wipe it dry immediately.

After cleaning and disinfections take care of the wooden parts with a special wood preservative.

10 Warranty

The warranty only concerns manufacturing defects or materials. It applies to the following items :

Bed frame	5 years
Electrical components	2 years (not including backup battery)

The warranty takes effect on the day on which the bed is purchased. If a component proves to be defective during the afore-mentioned period of validity, Thuisne shall undertake to replace or repair the said component unless the defect is due to incorrect handling.

Products that have been used incorrectly, neglected, poorly serviced or stored under inadequate conditions are not covered by the warranty. Although not covered is damage following to overload or violence.

Your bed has been supplied with adhesive stickers specifying the following serial numbers :

**Type plate (head part of the bed base and head part)
serial numbers on the back rest (head- and leg-raising devices)
and a serial number on each actuator**

The lifting pole has its own type plate.

Please quote these serial numbers in the event of a complaint.

11 Disposal of the various components at the end of the bed's useful life

Your Multicomfort II G2 bed consists of a welded, painted steel frame and electrical actuators. Thuisne urges you to comply with the waste disposal legislation currently in force in your country when your bed reaches the end of its useful life.

Do not dispose electrical components through the house garbage. Send them back to us. We dispose them properly.

12 What to do if the bed breaks down

You must always start by checking whether :

- the plug is inserted and there is tension in the network
- all the components are correctly connected to each other ;

- the blocking device is in the unlocked position (figure 23, pictogram lock opened) ;
- press firmly the green button on the plug once

If the problem persists, proceed as follows :

 Problem encountered		 Possible solution
What are the motors doing ?	Is there a switching noise in the plug when a function button is pressed	Measure
The motors are not functioning	no	<p>Measure: Press the green button on the mains plug. If the bed works again change the backup battery. Check if the locking switch is in open position.</p> <p>Other possible reasons</p> <ul style="list-style-type: none"> • Hand control defect • Electric faulty • Locking box defect. <p>Please contact your supplier!</p>
Only one motor is not functioning	yes	<p>• Defective motor or not correct connected to the main units</p> <p>Measure:</p> <ul style="list-style-type: none"> • Check the connection • Unplug and change the motor <p>Please contact your supplier!</p>
The motors have stopped or are not working	yes	<p>• Overload of the actuator due to high loads or to long operating time without pause.</p> <p>Measure:</p> <ul style="list-style-type: none"> • Unload the bed or stop operating for at least 10 minutes <p>For the bed base:</p> <p>Unload the bed base and if required bring it in the horizontal position (see 7.5)</p> <p>For the height adjustment:</p> <p>Never attempt to alter the height of the bed by removing a motor. Because the bed may be still operating by the backup battery for a short time when the overload security is working, we recommend to change the backup battery.</p>
One motor stops, although it has not reached its end of the stroke and the button of the hand control is pressed.	yes	<p>The actuator is overloaded or has been used for too long without a pause:</p> <p>Action:</p> <p>Stop operating for at least 10 minutes.</p> <p>Because the bed may be still operating by the backup battery for a short time when the overload security is working, we recommend to change the backup battery.</p>

If the problem resists, disconnect the drive system from the network and contact your specialised dealer, who delivered your care bed. Please describe him the encountered problem and inform him about the serial number of your bed.

Never open a motor yourself and do not use the bed if the electrical components are defective.


Without the expressed permission of Thuisne, the repair or the replacement of spare parts is allowed by a qualified technician only.

Any unauthorised intervention with the electrical components or any manipulation of the system will exclude the warranty. There may be danger for life by an electric shock because of improper intervention
We recommend in any case of a failure you have to contact your approved Thuisne retailer. The result of the failure must be stipulated so as to ensure an optimum solution to the problem.

Original – Please use copies only
Servicing schedule for the Thuisne care bed

Type :	Serial number :	Year of manufacture:
At least once a year after installation and re-installation (see local legislation, which may impose more controls)		
<u>Maintenance reason :</u>		
Annually	After repair	Before re-installation
Evaluation of the general electricity supply1 (mandatory in Germany and recommended by Thuisne for other countries) This check procedure must be carried out by a qualified electrician or trained electro-technician using appropriate measuring and testing equipment. Emergency fault current of the mains unit ≤0,100 mA (EN 60601-1) Insulation resistance ≥7 MOhm (see VDE 0751) Remark : the system is not working properly if the values deviate by more than 50% from those measured the first/previous time.		OK
Visual inspection :general condition		
Do all the moveable parts work easily ? Are there any parts that hang loosely or rattle when the bed is moved ? Are there any scraping or grinding noises when the bed is moved ? Are there any rusted parts ? Are the bed base, head and foot of the bed, backrest, seat, foot rest, bracket and castors twisted ? Are there any visible tears to the load-bearing components ? Are the motor supports intact or free from other visible signs of deterioration ? Have the side rails been altered ? Have they deteriorated or lost their shape ? Are the side rails easy to replace ? Do the castors turn easily and are the brakes working ? Is the plastic housing over the actuators damaged ? Is the transformer housing still firmly attached to the motor of the head-raising device ? Are all the covers intact, tear-free and firmly attached ? Are the plug and the inner compartment (cover plate at the upper end of the head) free from signs of corrosion and tears ? Is the lifting pole free from corrosion and tears ? Are the triangle and flange intact and tear-free ? Are all the loops present and intact ?		OK
Overall function check		
Do the side rails firmly engage into position when raised and can be released again if lifted up a little ? Do the side rails support a vertical force of 75 kg ? Once the side rails have been clicked into position, is there a maximum vertical gap of 10 mm ? Check that the deviation between the two tubes (Paris) and the two plates (Milano) has not been modified and is still less than 120 mm ?		OK
Actuator function check		
Drive through the complete the adjustment positions until the actuators completely stop and check the actuator range (variable height 400 to 800 mm, 110 mm for the head rest, 39 mm for the foot rest) Noise level, speed ? Does the stopping device work when all the actuators have reached the end of the range ? Maximum intensity value less than 1 A ? Is the backup-battery charged enough ? Do the actuators of the height adjustment reach the highest position/incline under max. load in less than 3 seconds distance ?		OK
Mains cable		
Is the main cable intact ? Check the function and condition of the anti-traction system on the mains cable (see figure 16) and the plug.		OK
All cables		
Check that the O-ring sealing rings on the plug are present and in good condition. Check that the plastic covers and plugs are intact (no flaws, crushing, tears, etc.) ; spiral sections not stretched. Check that the positioning and attachment of the cables is consistent with the cable installation technique (see chapter 6.5) ?		OK
Hand control unit and blocking device		
All buttons are moving smoothly and do not stick ? System is blocked in locked position ? System does not work in position I and II		OK
Documentation		
Are the instructions for use available ? Are all the adhesive labels firmly attached and legible ? Furthermore, for safety reasons, the following measurements must be recorded every five years : • Change of the lifting pool ; • Careful check of all actuators (no evidence of cable and connection deterioration, noise level). Change of the backup battery every 2 years.		OK
Description of defects; remarks:		
Checked by :		
Checked on:		
Checked at:		


Plattform MULTICOMFORT II G2

Pos.Nr.	Description	Art.code
1	Drive for backrest without transformator (IPX4)	8300210
	Transformer (IPX4) for system with integrated locking box	8300222
	Transformer (IPX4) for system with integrated locking box, United Kingdom version	8300226
No pic.	Pin for fast unlocking of back rest actuator	8300800
	Hand control with integrated locking device	8300276
No pic.	Key for locking device	8300273
2	Leg rest actuator (IPX4)	8300231
No pic.	Cable holder (plastics)	8400690
9 7	Frame for back rest, complete	8400680
3 5 6	Frame for leg rest, complete, including Rest for calves Rest for thighs	8400681
10	Raster für Fowler-Position, komplett	8400350
11	Gliders for side rails, complete, including pins for the gliders	8400701
12	Side rails Milano (wood, pair)	214927122
No pic.	Side rails Paris (metall, pair)	2'
13	Thumb screw	8

Stand 10/05

Head and foot part MULTICOMFORT II G2

Pos. Nr	Description	Model	Art.-Nr.
1	Wooden blend head-/footpart	Mantovani	8400620
	Wooden blend heas part	Haydn	8400660
	Wooden blend foot part	Haydn	8400650
No pic.	Locking pins for Milano side rail	all	8400700
2	Actuator head- or foot part	all	8300120
3 4	Heave pile head part	all	8400401
	Heave pile foot part	all	8400411
5	Heave pile top 70x40 cm	all	8300820
6	Wheel (100 mm)	all	8400500
No pic.	Triangle grip and belt	all	8200901
No pic.	Transport kit	all	8400100
No pic.	Manual	all	9019'
No pic.	Colour spray RAL1019	all	82007

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